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BD Vacutainer® Sodium Fluoride and Potassium Oxalate Plus Tube with BD Hemogard™ Safety Closure

BD Diagnostics
Preenalytical Systems

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Standards:	ISO 6710, EN14820
Linee Guida:	Clinical and Laboratory Standards Institute (CLSI; Formerly NCCLS): Tubes and Additives for Venous Blood Specimen Collection; Approved Guideline (6th Edition). Document GP-A6. Wayne, PA, USA, 2010.

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
Istruzioni d'uso

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- Guder WG, Narayanan S, Wisser H and Zawta B. Samples: From the Patient to the Laboratory: the Impact of Preanalytical Variables on the Quality of Laboratory Results (4th Edition). Darmstadt, Germany: Wiley-VCH; 2009.
- Evans MJ, Livesey JH, Ellis MJ and Yandle TG. "Effect of Anticoagulants and Storage Temperatures on Stability of Plasma and Serum Hormones". Clin Biochem. 2001; 34: 107-112.
- Chan AY, Swaminathan R and Cockram CS. "Effectiveness of Sodium Fluoride as a Preservative of Glucose in Blood". Clin Chem. Feb 1989; 35: 315 – 317.
- Uchida K, Matuse R, Toyoda E, Okuda S and Tomita S. "A New Method of Inhibiting Glycolysis in Blood Samples". Clinica Chimica Acta. Feb. 1988; 172(1): 101-8.0
- BD White Paper VS7595-OUS: "An Evaluation of Lactate Stability over 24 Hours at Room Temperature in Various BD Vacutainer® Blood Collection Tubes Containing Glycolytic Inhibitor/Anticoagulant Combinations". 2008.
- BD White Paper VS5918: "BD Vacutainer® PLUS Sodium Fluoride/Potassium Oxalate Tubes: Evaluation of Clotting, Hemolysis and Glucose at Initial Time of Collection and After 4 and 24 Hours Storage at Room Temperature When Compared to Glass Sodium Fluoride Tubes". 2008.

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Bibliografia

1. Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS): Procedures for the Handling and Processing of Blood Specimens; Approved Guideline (4th Edition). Document H18-A4. Wayne, PA, USA: 2010.
2. Guder WG, et al. Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine for Quality of Diagnostic Samples (3rd Edition). Darmstadt, Germany: GIT, 2010.
3. Tietz NW. Clinical Guide to Laboratory Tests (4th Edition). W.B. Saunders, USA: 2006.

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1. Guder WG, Narayanan S, Wisser H and Zawta B. Samples: From the Patient to the Laboratory: the Impact of Preanalytical Variables on the Quality of Laboratory Results (4th Edition). Darmstadt, Germany: Wiley-VCH; 2009.
2. Evans MJ, Livesey JH, Ellis MJ and Yandle TG. "Effect of Anticoagulants and Storage Temperatures on Stability of Plasma and Serum Hormones". Clin Biochem. 2001; 34: 107-112.
3. Chan AY, Swaminathan R and Cockram CS. "Effectiveness of Sodium Fluoride as a Preservative of Glucose in Blood". Clin Chem. Feb 1989; 35: 315 - 317.
4. Uchida K, Matusz R, Toyoda E, Okuda S and Tomita S. "A New Method of Inhibiting Glycolysis in Blood Samples". Clinica Chimica Acta. Feb. 1988; 172(1): 101-8.0
5. BD White Paper VS7595-0US: "An Evaluation of Lactate Stability over 24 Hours at Room Temperature in Various BD Vacutainer® Blood Collection Tubes Containing Glycolytic Inhibitor/Anticoagulant Combinations". 2008.
6. BD White Paper VS7594: "A Comparative Evaluation of BD Vacutainer® Plus Fluoride/EDTA Tubes for Glucose & HbA1c Preservation After 24 Hours Storage at Room Temperature". 2007.
7. BD White Paper VS5905: "BD Vacutainer® PLUS Sodium Fluoride/EDTA Tubes: Evaluation of Glucose and Hemoglobin A1C at Initial Time of Collection and After 24 Hour Storage at 4°C". 2001.
8. BD White Paper VS5904: "BD Vacutainer® Plus Sodium Fluoride EDTA: Evaluation of Clotting, Haemolysis and Glucose at Initial Time of Collection and After 24h Storage at RT When Compared to Glass Sodium Fluoride/Potassium Oxalate Tubes". 2001

BD Diagnostics
Preanalytical Systems
Becton Dickinson Italia SpA
Via Enrico Cialdini, 16
20161 MILANO (MI)
Tel: 02.48240.1
Fax: 02.48240344
www.bd.com
preanalytical@bd.europe.com

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BD Vacutainer® Sodium Fluoride and Sodium EDTA Plus Tube with BD Hemogard™ Safety Closure

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TIS367933, AC, CF, 23 05 14, 01

BD Diagnostics
Preanalytical Systems

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- Guder WG, Narayanan S, Wisser H and Zawta B. Samples: From the Patient to the Laboratory: the Impact of Preanalytical Variables on the Quality of Laboratory Results (4th Edition). Darmstadt, Germany: Wiley-VCH; 2009.
- Evans MJ, Livesey JH, Ellis MJ and Yandle TG. "Effect of Anticoagulants and Storage Temperatures on Stability of Plasma and Serum Hormones". Clin Biochem. 2001; 34: 107-112.
- Chan AY, Swaminathan R and Cockram CS. "Effectiveness of Sodium Fluoride as a Preservative of Glucose in Blood". Clin Chem. Feb 1989; 35: 315 - 317.
- Uchida K, Matusz R, Toyoda E, Okuda S and Tomita S. "A New Method of Inhibiting Glycolysis in Blood Samples". Clinica Chimica Acta. Feb. 1988; 172(1): 101-8.0
- BD White Paper VS7595-OUS: "An Evaluation of Lactate Stability over 24 Hours at Room Temperature in Various BD Vacutainer® Blood Collection Tubes Containing Glycolytic Inhibitor/Anticoagulant Combinations". 2008.
- BD White Paper VS7594: "A Comparative Evaluation of BD Vacutainer® Plus Fluoride/EDTA Tubes for Glucose & HbA1c Preservation After 24 Hours Storage at Room Temperature". 2007.
- BD White Paper VS5905: "BD Vacutainer® PLUS Sodium Fluoride/EDTA Tubes: Evaluation of Glucose and Hemoglobin A1C at Initial Time of Collection and After 24 Hour Storage at 4°C". 2001.
- BD White Paper VS5904: "BD Vacutainer® Plus Sodium Fluoride EDTA: Evaluation of Clotting, Haemolysis and Glucose at Initial Time of Collection and After 24h Storage at RT When Compared to Glass Sodium Fluoride/Potassium Oxalate Tubes". 2001

Bibliografia

- Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS): Procedures for the Handling and Processing of Blood Specimens; Approved Guideline (4th Edition). Document H18-A4. Wayne, PA, USA: 2010.
- Guder WG, et al. Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine for Quality of Diagnostic Samples (3rd Edition). Darmstadt, Germany: GIT, 2010.
- Tietz NW. Clinical Guide to Laboratory Tests (4th Edition). W.B. Saunders, USA: 2006.

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