Nova Pro Glucose Test Strips - For Use ONLY with the Nova Pro Family of Meters

Nova Pro Teststreifen zur ausschließlichen Verwendung mit Messgeräten der Nova Pro-Familie, Tiras reactivas para glucosa Nova Pro - Para uso EXCLUSIVO con la familia de medidores Nova Pro, Bandelettes réactives de glucose Nova Pro - À utiliser UNIQUEMENT avec la famille de lecteurs de glycémie Nova Pro, Strisce reattive al glucosio Pro Nova - Per uso ESCLUSIVO con la Serie di misuratori Pro Nova, Tiras de análise da glicose Nova Pro - Para utilização APENAS com a Família de Medidores Nova Pro

Table 1 Typical Within-run Precision - Blood Results, Tabelle 1 Typische Blutproben-Messgenauigkeit, Tabla 1 Precisión intraserial típica – Resultados de análisis de sangre, Tableau 1 Précision typique de l'analyse – Résultats de sang, Tabella 1 Precisione tipica intra-analisi – Risultati con un campione di sangue, Tabela 1 Precisão típica dentro dos valores previsíveis – Resultados da análise ao sangue

| <div><div></div><div>Glu</div></div> LEVEL 1 | n=60 | <div><div></div><div>Glu</div></div> LEVEL 2 | n=60 | <div><div></div><div>Glu</div></div> LEVEL 3 | n=60 | <div><div></div><div>Glu</div></div> LEVEL 4 | n=60 | <div><div></div><div>Glu</div></div> LEVEL 5 | n=60 |
|--|------|--|------|--|------|--|------|--|------|
| <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | |
| 69.4 mg/dL 2.7 | | 155.7 mg/dL 2.1 | | 232.5 mg/dL 2.7 | | 337.1 mg/dL 2.9 | | 545.6 mg/dL 3.6 | |
| 8.6 mmol/L | | 8.6 mmol/L | | 12.9 mmol/L | | 18.7 mmol/L | | 30.3 mmol/L | |

Table 2 Typical Within-run Precision - Linearity Solution Results, Tabelle 2: Typische Linearitätslösung-Messgenauigkeit, Tabla 2 Precisión intraserial típica – Resultados de soluciones de linealidad, Tableau 2 Précision typique de l'analyse – Résultats de solution de linéarité, Tabella 2 Precisione tipica intra-analisi – Risultati con la soluzione di linearità, Tabela 2 Precisão típica dentro dos valores previsíveis – Resultados da solução de linearidade

| <div><div></div><div>LINEARITY</div></div> LEVEL 1 | n=60 | <div><div></div><div>LINEARITY</div></div> LEVEL 2 | n=60 | <div><div></div><div>LINEARITY</div></div> LEVEL 3 | n=60 | <div><div></div><div>LINEARITY</div></div> LEVEL 4 | n=60 | <div><div></div><div>LINEARITY</div></div> LEVEL 5 | n=60 |
|--|------|--|------|--|------|--|------|--|------|
| <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | | <div><div></div><div>X</div></div> CV(%) | |
| 20.9 mg/dL 4.4 | | 64.4 mg/dL 2.5 | | 103.1 mg/dL 2.7 | | 299.5 mg/dL 5.3 | | 484 mg/dL 1.9 | |
| 1.2 mmol/L 4.4 | | 3.6 mmol/L 2.5 | | 5.7 mmol/L 2.7 | | 16.6 mmol/L 5.3 | | 26.9 mmol/L 1.9 | |

EN

Contents

53881: Contains 100 strips (50 strips x 2 vials).

Intended Use

The Nova Pro Glucose Test Strip is intended for use only with the Nova Pro Family of Meters to quantitatively measure glucose (sugar) in whole blood. The Nova Pro Glucose Meter is plasma calibrated to allow easy comparison of results with laboratory methods. Nova Pro Glucose Test Strips are for testing outside the body (*in vitro* diagnostic use). It is not for diagnosis of or screening for diabetes.

Summary and Principle of Test

The Nova Pro Glucose Test Strip is designed with an electrode that measures glucose levels. Glucose (sugar) in the blood sample mixes with reagent on the test strip that Produces an electric current. The amount of current that is Produced depends on how much glucose is in the blood. Nova Pro Glucose Meter with Nova Pro Glucose Test Strip measures the glucose level in blood. The glucose result is displayed on the screen of the meter.

Glucose Interferences:

The Nova Pro Glucose Meter exhibits no interference from the following substances at known therapeutic levels: Acetaminophen, Ascorbic Acid, Dopamine, Ephedrine, N-acetylcysteine, D(+) Galactose, Ibuprofen, L-Dopa, Methyl-Dopa, Salicylate, Tetracycline, Tolazamide, and Tobutamide.

The Nova Pro Glucose and β -Ketone Meter exhibits no interference from the following substances at or above the upper clinical normal range concentrations: Bilirubin, Cholesterol, Creatinine, Triglycerides, and Uric Acid.

The Nova Pro Glucose Meter exhibits no interference from the following substances at the normal therapeutic levels found in renal dialysis: D Maltese

The Nova Pro Glucose Meter exhibits no interference in blood specimens with hematocrits from 20% to 65% or with varying oxygen content.

Reagent Composition

Each Nova Pro Glucose Test Strip contains the enzyme glucose dehydrogenase GDH-FAD (Aspergillus sp.) greater than or equal to 0.6 IU; additional ingredients (mediator, buffer, etc.) greater than or equal to 0.16 μ g. The test strip vial contains a molecular sieve desiccant.

Precautions

- Prior to use, read the Nova Pro Glucose Meter Instructions for Use Manual.
- Capillary blood glucose testing may not be appropriate for persons with decreased peripheral blood flow, as it may not reflect the true physiological state. Examples include, but are not limited to, severe hypotension, shock, hyperosmolar-hyperglycemia (with or without ketosis) and severe dehydration.
- DO NOT reuse test strips. Strips are designed to be disposed after a single use.
- Discard used test strips according to local regulations.
- Use only the Nova Pro Glucose Test Strips for testing with the Nova Pro Glucose Meter.
- If test result is higher or lower than expected, run a control solution test.
- Remove the test strip from the vial only when ready to test.
- Do not use the test strip if the expiration date has passed, for this may cause inaccurate results.
- Do not tamper with the test strip.

Storage and Handling

Keep the Nova Pro Glucose Test Strip vial tightly closed when not in use. Test strips should be stored only in original vial. Store the test strips at a temperature of 1-30°C. Do not freeze the test strips.

Expiration

The expiration date is printed on the vial of test strips. Once opened, the Nova Pro Test Strips are stable when stored as indicated for up to 180 days or until the expiration date, whichever comes first.

Test Procedure

Refer to your Nova Pro Glucose and β -Ketone Meter Instructions for Use Manual for detailed meter information prior to testing. Items Needed: Nova Pro Glucose and β -Ketone Meter, lancet, and Nova Pro Glucose Test Strips.

Normal Glucose Values:

The normal adult fasting blood glucose range for a person is 70-110 mg/dL (3.9-6.1 mmol/L).¹

High Glucose Values:

If the test result is higher than 600 mg/dL (33.3 mmol/L), the Glucose Meter will read "HI." Retest using a new test strip.

Low Glucose Values:

If the test result is lower than 10 mg/dL (0.6 mmol/L), the Glucose Meter will read "LO." Retest using a new test strip.

If the same HI or LO reading is obtained, run a control test.

Screen Alerts:

The Nova Pro Glucose Meter displays screen alerts.

- Battery Low - Change the battery or place the meter onto the Charging Station.
- Temperature Error - Meter will only work within the temperature range of 5°C to 40°C.
- Strip Removed - The test strip was removed before completing the test. The test is canceled.
- Wrong Strip - Insert the correct strip.
- Short Sample - Not enough sample was put onto the strip. Redo the test with a new strip.

Quality Control: How to Check the System

Run Control Solution Test:

The Nova Pro Glucose and β -Ketone Control Solution is used as a quality control check to make sure that the Nova Pro Glucose and β -Ketone Meter and the Nova Pro Test Strips are working correctly. The Nova Pro Control Solution is intended for use only with the Nova Pro Glucose and β -Ketone Test Strips.

The control solution contains a measured amount of glucose that reacts with the Nova Pro Glucose Test Strips. This reaction is similar to the way glucose in the blood reacts with the test strip. A control solution test is performed the same way that a blood glucose test is performed except that control solution is used in place of a blood drop.

The control solution test confirms that the meter and test strips are working correctly. The control solution test results should fall within the range of results printed on the vial label of the control.

When to Perform a Quality Control Test

Run 2 different levels of control during each 24 hours of testing and under the following conditions:

- As required by your institution's quality control policy or local regulatory requirements.
- If a patient test has been repeated and the blood glucose results are still lower or higher than expected.
- If you leave the test strip vial open for an extended period of time.
- If you drop the meter.

The Nova Pro Glucose and β -Ketone Control Solution should produce results that fall within the range of results printed on the control solution vial label of the control being used. If the control solution test result is outside the range (is either higher or lower), the Nova Pro Glucose Meter and test strip may not be working as a system.

Reasons for Out-of-Range Control Result:

- You may not be doing the test Properly. Retest and follow the instructions carefully.
- The control solution may have expired or have been contaminated. Check the expiration date on the control solution vial. Control solution is good for only 3 months after opening. Make sure the control solution vial is closed when not in use.
- The test strip may have expired. Check the expiration date on the test strip vial.
- The test strip may have been damaged. Retest using a new test strip.
- The Nova Pro Glucose and β -Ketone Meter may not be working Properly.

Limitations

The Nova Pro Glucose and β -Ketone Meter should have results that agree with a laboratory result to within 20%.² A test result within this range is considered accurate when testing with the Nova Pro Meter.

See the Performance Characteristics section for accuracy and precision information. However, there are factors that may cause results to differ by more than 20% in some situations. These factors are listed below:

- Blood source - Use only whole blood. Do not use serum or plasma.
- Venous and capillary blood may differ in glucose concentration by as much as 70 mg/dL (3.9 mmol/L), depending on the time of blood collection after food intake.¹
- Temperature and humidity extremes - Test results are best obtained when Nova Pro Systems are used within an operating relative humidity of 10-90% (non-condensing). Testing outside these ranges may cause inaccurate results.
- Altitude - There is no effect of altitudes up to 11,000 feet (3353 meters) above sea level.
- Venous Specimens - Fluoride (gray top tubes) or EDTA (purple top tubes) should not be used as preservatives for venous specimens.

Additional Information for Healthcare Professionals

- The Nova Pro Glucose Test Strips are calibrated against plasma.
- The Nova Pro Glucose Test Strips are used for whole blood testing only.
- Healthcare Professionals and others using this system on multiple patients should be aware that all Products or objects that come into contact with human blood should be handled as if capable of transmitting viral diseases, even after cleaning.

Performance Characteristics

The test result range for the Nova Pro Glucose and β -Ketone Meter is 10 to 600 mg/dL (0.6 to 33.3 mmol/L) for Glucose.

Accuracy-Nova Pro Meter

The accuracy of the Nova Pro Glucose and β -Ketone Meter was assessed by comparing the glucose results obtained from the Nova Pro System with the results obtained from the YSI 2300 Stat Plus Analyzer (Yellow Springs Instrument Co., Yellow Springs, OH) and the Nova StatStrip Xpress Blood Glucose Monitoring System. Specimens were obtained and tested at 3 sites within one hospital.

| | Accuracy | | |
|-----------------------------|-----------|--------|----------|
| | Capillary | Venous | Arterial |
| Number of Samples | 153 | 125 | 125 |
| Slope | 0.960 | 0.987 | 0.972 |
| Y-intercept | -0.067 | -1.778 | 2.907 |
| R-squared (r ²) | 0.9961 | 0.9944 | 0.9935 |
| Compared to | YSI | YSI | YSI |

Typical Precision

Precision of the Nova Pro Glucose Meter and β -Ketone System was measured with both spiked whole blood and control solutions in the laboratory. The typical strip lot results are shown in Table 1 and Table 2.

REFERENCES

- Burtis, Carl A. and Ashwood, Edward R., ed. 1998. Carbohydrates (Sachs, DB): Tietz Textbook of Clinical Chemistry. Philadelphia, PA: W. B. Saunders Co.
- Clarke, W.L., Cox, D., Gonder-Frederich, L.A., et al. 1987. Evaluating Clinical Accuracy of Systems for Self monitoring of Blood Glucose. Diabetes Care. Vol 10 pp 622-628.

CONTACT IN THE US: Nova Biomedical Technical Support at 781.894.0800 or toll-free 1.800.545.6682.

CONTACT IN CANADA: 1.800.263.5999.

INTERNATIONAL: Contact your Nova Biomedical Distributor located on the web site www.novabiomedical.com.

NOT FOR EMERGENCY OR MEDICAL INFORMATION
Labeling Symbols Used: Refer to the Instructions for Use Manual.

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