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Bibliografia

1. Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS): Procedures for the Handling and Processing of Blood Specimens; Approved Guideline (4th Edition). Document H18-A4. Wayne, PA, USA: 2010.
2. Guder WG, et al. Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine for Quality of Diagnostic Samples (3rd Edition). Darmstadt, Germany: GIT, 2010.
3. Tietz NW. Clinical Guide to Laboratory Tests (4th Edition). W.B. Saunders, USA: 2006.

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1. Giavarina D, Fortunato A, Barzon E, Church S, Bérubé J, Green S and Siffiati G. "Evaluation of BD Vacutainer® PST™ II Tubes for a Wide Range of Immunoassays". Clin Chem Lab Med. 2009; 47(2): 237-41.
2. Chance J, Bérubé J, Vandersmissen M, Blanckaert N. "Evaluation of the BD Vacutainer® PST™ II Tube for Special Chemistry Analytes". Clin Chem Lab Med, 2009; 47(3): 358-61.
3. van den Ouweland JMW, Church S. "High Total Protein Impairs Appropriate Gel Barrier Formation in BD Vacutainer Blood Collection Tubes". Clin Chem. 2007; 53(2): 364-5.
4. Morovat A, James TS, Cox SD, Norris SG, Rees MC, Gales MA and Taylor RP. "Comparison of Bayer Advia Centaur® Immunoassay Results Obtained on Samples Collected in Four Different Becton Dickinson Vacutainer® Tubes". Ann Clin Biochem. 2006; 43: 481-7.
5. Srivasta R, Murphy MJ, Card J, Severn A and Fraser C.G. "The Case of the Floating Gel". J Clin Pathol. 2004; 57: 1333-4.
6. Dimeski G, Badrick T, Flatman R and Ormiston B. "Roche IFCC Methods for Lactate Dehydrogenase Tested for Duplicate Errors with Greiner and Becton-Dickinson Lithium-Heparin and Greiner Serum Samples". Clin Chem. 2004; 5-(12): 2391-92.
7. Spiritus T, Zaman Z and Desmet W. "Iodinated Contrast Media Interfere with Gel Barrier Formation in Plasma and Serum Separator Tubes". Clin Chem. 2003; 49(7): 1187-9.
8. Lin F, Cohen R, Losada R and Bush V. "Cellular Sedimentation and Barrier Formation under Centrifugal Force in Blood Collection Tubes". Lab Med. Oct. 2001; 32(10).
9. BD White Paper VS7597-OUS: "A Comparative Evaluation of BD Vacutainer® PST™ II Tubes with BD Vacutainer® Lithium Heparin Plus and BD Vacutainer® Serum Plus Tubes for Selected Hormones, Therapeutic Drugs, Tumor Markers and Other Chemistry Analytes". 2008.
10. BD White Paper VS7513-OUS: "Performance of BD Vacutainer® PST™ II Tubes at Four and Five Minute Centrifugation Times". 2006.
11. BD White Paper VS7046: "Performance of BD Vacutainer® PST™ II PLUS Tubes at Recommended Centrifugation Conditions". 2002.
12. BD White Paper VS5925: Analyte Stability Supports Extended Use of Plasma Collected in BD Vacutainer® PST™ II Plastic Tubes. 2001.

BD Diagnostics
Preanalytical Systems
Becton Dickinson Italia SpA
Via Enrico Cialdini, 16
20161 MILANO (MI)
Tel: 02.48240.1
Fax: 02.48240344
www.bd.com
preanalytical@bd.europe.com

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BD Vacutainer® PST™ II Plus Tube with BD Hemogard™ Safety Closure

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TIS368497, CF, 23 05 14, 01

BD Diagnostics
Preanalytical Systems

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Standards: ISO6710, EN14820

Linee Guida: Clinical and Laboratory Standards Institute (CLSI; Formerly NCCLS):
Tubes and Additives for Venous Blood Specimen Collection; Approved
Guideline (6th Edition). Document GP-A6. Wayne, PA, USA, 2010.

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
1. Giavarina D, Fortunato A, Barzon E, Church S, Bérubé J, Green S and Siffati G. "Evaluation of BD Vacutainer® PST™ II Tubes for a Wide Range of Immunoassays". Clin Chem Lab Med. 2009; 47(2): 237-41.
2. Chance J, Bérubé J, Vandersmissen M, Blanckaert N. "Evaluation of the BD Vacutainer® PST™ II Tube for Special Chemistry Analytes". Clin Chem Lab Med, 2009; 47(3): 358-61.
3. van den Ouweland JMW, Church S. "High Total Protein Impairs Appropriate Gel Barrier Formation in BD Vacutainer Blood Collection Tubes". Clin Chem. 2007; 53(2): 364-5.
4. Morovat A, James TS, Cox SD, Norris SG, Rees MC, Gales MA and Taylor RP. "Comparison of Bayer Advia Centaur® Immunoassay Results Obtained on Samples Collected in Four Different Becton Dickinson Vacutainer® Tubes". Ann Clin Biochem. 2006; 43: 481-7.
5. Srivasta R, Murphy MJ, Card J, Severn A and Fraser C.G. "The Case of the Floating Gel". J Clin Pathol. 2004; 57: 1333-4.
6. Dimeski G, Badrick T, Flatman R and Ormiston B. "Roche IFCC Methods for Lactate Dehydrogenase Tested for Duplicate Errors with Greiner and Becton-Dickinson Lithium-Heparin and Greiner Serum Samples". Clin Chem. 2004; 5-(12): 2391-92.
7. Spiritus T, Zaman Z and Desmet W. "Iodinated Contrast Media Interfere with Gel Barrier Formation in Plasma and Serum Separator Tubes". Clin Chem. 2003; 49(7): 1187-9.
8. Lin F, Cohen R, Losada R and Bush V. "Cellular Sedimentation and Barrier Formation under Centrifugal Force in Blood Collection Tubes". Lab Med. Oct. 2001; 32(10).
9. BD White Paper VS7597-OUS: "A Comparative Evaluation of BD Vacutainer® PST™ II Tubes with BD Vacutainer® Lithium Heparin Plus and BD Vacutainer® Serum Plus Tubes for Selected Hormones, Therapeutic Drugs, Tumor Markers and Other Chemistry Analytes". 2008.
10. BD White Paper VS7513-OUS: "Performance of BD Vacutainer® PST™ II Tubes at Four and Five Minute Centrifugation Times". 2006.
11. BD White Paper VS7046: "Performance of BD Vacutainer® PST™ II PLUS Tubes at Recommended Centrifugation Conditions". 2002.
12. BD White Paper VS5925: Analyte Stability Supports Extended Use of Plasma Collected in BD Vacutainer® PST™II Plastic Tubes. 2001.

BD Diagnostics
Preanalytical Systems
Becton Dickinson Italia SpA
Via Enrico Cialdini, 16
20161 MILANO (MI)
Tel: 02.48240.1
Fax: 02.48240344
www.bd.com
preanalytical@bd.europe.com