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## BD Vacutainer® Lithium Heparin Plasma Plus Tube with BD Hemogard™ Safety Closure

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TIS368886, CF, 23 05 14, 02

BD Diagnostics  
Preanalytical Systems

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Standards:	ISO6710, EN14820
Linee Guida:	Clinical and Laboratory Standards Institute (CLSI; Formerly NCCLS): Tubes and Additives for Venous Blood Specimen Collection; Approved Guideline (6th Edition). Document GP-A6. Wayne, PA, USA, 2010.

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
## Istruzioni d'uso

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- Guder WG, Narayanan S, Wissner H and Zawta B. Samples: From the Patient to the Laboratory: the Impact of Preanalytical Variables on the Quality of Laboratory Results (4th Edition). Darmstadt, Germany: Wiley-VCH; 2009.
- Giavarina D, Fortunato A, Barzon E, Church S, Bérubé J, Green S and Siffiat G. "Evaluation of BD Vacutainer® PST™ II Tubes for a Wide Range of Immunoassays". Clin Chem Lab Med. 2009; 47(2): 237-41.
- Chance J, Bérubé J, Vandersmissen M, Blanckaert N. "Evaluation of the BD Vacutainer® PST™ II Tube for Special Chemistry Analytes". Clin Chem Lab Med, 2009; 47(3): 358-61.
- Nosanchuk JS, Stull R and Keefner R. "The Effects of Substitution of Plasma for Serum on Chemistry Stat Turnaround Time". Lab Med. July 1991; 22(7): 465-9.0
- Harr R, Bond L and Trumbull D. "A Comparison of Results for Serum versus Heparinized Plasma for 30 Common Analytes". Lab Med. July 1987; 18(7): 449-55.
- BD White Paper VS7597-0US: "A Comparative Evaluation of BD Vacutainer® PST™ II Tubes with BD Vacutainer® Lithium Heparin Plus and BD Vacutainer® Serum Plus Tubes for Selected Hormones, Therapeutic Drugs, Tumor Markers and Other Chemistry An
- BD White Paper VS7174: "Comparison of BD Vacutainer® Lithium Heparin Plus Tubes with BD Vacutainer® Lithium Heparin Glass Tubes for Routine Chemistry Analytes". 2004.

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## Bibliografia

- Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS): Procedures for the Handling and Processing of Blood Specimens; Approved Guideline (4th Edition). Document H18-A4. Wayne, PA, USA: 2010.
- Guder WG, et al. Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine for Quality of Diagnostic Samples (3rd Edition). Darmstadt, Germany: GiT, 2010.
- Tietz NW. Clinical Guide to Laboratory Tests (4th Edition). W.B. Saunders, USA: 2006.

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Standards:	ISO 6710, EN14820
Linee Guida:	Clinical and Laboratory Standards Institute (CLSI; Formerly NCCLS): Tubes and Additives for Venous Blood Specimen Collection; Approved Guideline (6th Edition). Document GP-A6. Wayne, PA, USA, 2010.

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
<sup>2,3</sup>

## Bibliografia

1. Clinical and Laboratory Standards Institute (CLSI; formerly NCCLS): Procedures for the Handling and Processing of Blood Specimens; Approved Guideline (4th Edition). Document H18-A4. Wayne, PA, USA: 2010.
2. Guder WG, et al. Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine for Quality of Diagnostic Samples (3rd Edition). Darmstadt, Germany: GiT, 2010.
3. Tietz NW. Clinical Guide to Laboratory Tests (4th Edition). W.B. Saunders, USA: 2006.

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1. Guder WG, Narayanan S, Wissner H and Zawta B. Samples: From the Patient to the Laboratory: the Impact of Preanalytical Variables on the Quality of Laboratory Results (4th Edition). Darmstadt, Germany: Wiley-VCH; 2009.
2. Giavarina D, Fortunato A, Barzon E, Church S, Bérubé J, Green S and Siffiati G. "Evaluation of BD Vacutainer® PST™ II Tubes for a Wide Range of Immunoassays". Clin Chem Lab Med. 2009; 47(2): 237-41.
3. Chance J, Bérubé J, Vandersmissen M, Blanckaert N. "Evaluation of the BD Vacutainer® PST™ II Tube for Special Chemistry Analytes". Clin Chem Lab Med, 2009; 47(3): 358-61.
4. Nosanchuk JS, Stull R and Keefner R. "The Effects of Substitution of Plasma for Serum on Chemistry Stat Turnaround Time". Lab Med. July 1991; 22(7): 465-9.0
5. Harr R, Bond L and Trumbull D. "A Comparison of Results for Serum versus Heparinized Plasma for 30 Common Analytes". Lab Med. July 1987; 18(7): 449-55.
6. BD White Paper VS7597-0US: "A Comparative Evaluation of BD Vacutainer® PST™ II Tubes with BD Vacutainer® Lithium Heparin Plus and BD Vacutainer® Serum Plus Tubes for Selected Hormones, Therapeutic Drugs, Tumor Markers and Other Chemistry An
7. BD White Paper VS7174: "Comparison of BD Vacutainer® Lithium Heparin Plus Tubes with BD Vacutainer® Lithium Heparin Glass Tubes for Routine Chemistry Analytes". 2004.

**BD Diagnostics**  
Preanalytical Systems  
Becton Dickinson Italia SpA  
Via Enrico Cialdini, 16  
20161 MILANO (MI)  
Tel: 02.48240.1  
Fax: 02.48240344  
www.bd.com  
preanalytical@bd.europe.com