cryokit

criogenic protection

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Storage / Storage, Use, Cleaning, Maintenance, and Disinfection Instructions

The apron and gaiters, both new and used, must be examined in advance, before each use, in order to ensure that there are no conditions of wear, abrasions, lacerations or cuts, such as to compromise the degree of protection required. Any contaminants must be immediately removed. Any contamination, alteration or incorrect use of the device dangerously reduces its performance. Store in a cool dry place away from light.

The results of the tests on the prototypes determined a level / class of protection provided by the PPE

PPE class 3

Accessories that can be used together with PPE

Gloves, aprons and visors

Protection classes referring to the different risk levels and the corresponding limits of use.

The apron and leggings do not protect against:

- Chemical agents other than those indicated
- Ionizing radiations
- Biological agents
- Any risk other than that indicated in this information note

Duration and / or expiry of the PPE / component (month / year), if applicable

Each apron or gaiter is contained in a polythene envelope with this information note inside. To waste of the envelope refer to the provisions of the municipality to which it belongs.

Risks against which the PPE is designed

Processes that require protection against splashes of liquid nitrogen and cryogenic gases in general in clinical laboratories, the pharmaceutical industry, cryopreservation, bio-medical, biological banks, the food industry, handling of dry ice, etc. Possible leakage of liquid in the transfer operations (splash).

Name, address and identification number of the notified body in charge of assessing the PPE (release of Forms B, C2, D) *

ANCCP CERTIFICATION AGENCY srl - Via dello Struggino, 6
I-57121 Livorno (Italy) Tel: +39.0586.209006 Fax: +39.0586.278450
e-mail: anccp@anccp.it - n° Notified Body 0302
Harmonized standards / techniques used
EN388: 2016 / EN511: 2006

Website (where the EU Declaration of Conformity can be accessed)
www.cryokit.net