NOTA INFORMATIVA cod. 348085 - mod. Neo Eco

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INFORMATIVE NOTE cod. 348085 - mod. Neo Eco

Read this briefing note through carefully prior to use and prior to any maintenance work. The purpose of the information contained herein is to assist and direct the user in choosing and using the PPE. Neither the manufacturer nor the distributor can be held liable for an incorrect use of the PPE. This briefing note must be kept for the duration of the PPE. Glove model: **cod. 348085 - mod. Neo Eco**.

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The CE 0321 mark means that these gloves are personal protective equipment conforming to the essential health and safety requirements set forth in the (EU) Regulation 2016/425 and have been certified by the notified body: **SATRA Technology Centre Ltd., Wyndham Way, Telford Way, Kettering, Northamptonshire, NN16 8SD, United Kingdom, n° 0321**. In addition, the manufacturer is subject to the control of production CE quality guarantee by the notified body: **SATRA Technology Centre Ltd., Wyndham Way, Telford Way, Kettering, Northamptonshire, NN16 8SD, United Kingdom, n° 0321**.

SIZES

OUR SIZE REFERENCES	SIZES CONTEMPLATED EN 420:2003+AI:2009
M	8
L	9
XL	10

The dexterity level of the fingers is equivalent to 5 out of 5.  
The gloves meet the safety characteristics only if the size is right and worn properly.

RANGE OF USE

The gloves have been made as protection against mechanical risks such as abrasion, cutting, ripping and perforation and they also provide protection for the hands when coming into contact with chemicals and micro-organisms (intended as fungus, bacteria and viruses).

PICTOGRAMS AND PROTECTION LEVELS

EN 388: 2016	EN ISO 374-1:2016 - Type A	EN ISO 374-5:2016
 <b>2121X</b>	 <b>JKLNPT</b>	 <b>VIRUS</b>
<b>MECHANICAL RISKS</b>  MAX	<b>PROTECTION AGAINST CHEMICALS</b>	<b>PROTECTION AGAINST MICRO-ORGANISMS</b> Resistance to air and water penetration using the method referred to in EN 374-2:2014 standards
2 Abrasion 4 1 Cutting 5 2 Ripping 4 1 Perforation 4		
X: Cut resistance as per EN ISO 13997:1999. The glove has not been tested for this characteristic, as it is not applicable.		

CHEMICAL RESISTANCE TABLE according to EN ISO 374-1: 2016

Chemical	Permeation index (0-6)	Time of Passage (min)	Index of degradation (EN 374-4: 2013)
J N-heptane	2	30	78.5%
K Hydroxide of sodium 40%	6	480	-4.3%
L Sulphic acid 96%	4	120	34.9%
N Acetic acid 99%	3	60	34.9%
P Hydrogen peroxide 30%	6	480	0.3%
T Formaldehyde 37%	6	480	-1.0%

Time of passage. The time a chemical takes to pass through the glove material. The values may differ when the chemical substance is used in a compound. Do not use gloves in contact with chemicals any longer than the time specified above. Since the performance levels are based on lab test results, they may differ from the actual conditions at the workplace. Consequently it is advisable to conduct a preliminary test of the gloves to see if they are suitable for their intended use. The protective characteristics refer to new PPE, in good condition and that has never been subject to cleaning treatments.

USE

Prior to use, make sure the gloves are in good condition by which we mean no cuts, holes or cracks, etc. If the gloves are not in perfect condition, replace them immediately. The gloves must be used solely for the risks contemplated in this briefing note. Avoid using PPE near moving parts in which it could get entangled. Ensure that the chemical substances cannot get inside the gloves via the over sleeves. Movements, tears, frictions and degradation caused by contact with the chemical substance may significantly reduce the actual time of use. In the case of corrosive agents, degradation may be the most important factor to consider when choosing the gloves.

CLEANING

Do not leave the gloves in a contaminated condition if you intend using them again. It is advisable to wash them before taking them off to avoid all possible contamination of hands. Dry in the air. When the substances cannot be removed, take the gloves off alternatively, first the right then the left to avoid the substances touching bare hands and then dispose of the gloves.

STORING

The gloves are in polythene bags. Keep the gloves in their original packaging in a clean, dry place, away from sources of heat and direct sunlight. If stored as indicated, the gloves will keep their characteristics for a long time. The user must always inspect the gloves to ensure they are sound prior to using them. How long the gloves will last depends on the use and care the user has.

NOTE

Used gloves can be contaminated with chemicals or biological substances: dispose of them in compliance with the local regulations in force on this subject (dump, incinerator). The mixture the gloves are made with contains small quantities of rubber which could give rise to allergic reactions in people sensitive to this substance. In the case of manufacturing defects this PPE will be replaced.  
\* **In the event of discrepancies between the different translations, the Italian version will exclusively apply.**  
For more information, contact:

Registered Community Trade Mark, no 016928426 c/o EUIPO - Alicante - Spain

The EU conformity declaration is available at following address:  
[www.nerispa.com](http://www.nerispa.com)

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