Aesculap Dental Implantology

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Dentistry is currently facing drastic changes. The current state of medical and technological knowledge means that operations can be performed with ever greater precision and using less and less invasive techniques. Along with the availability of increasingly sophisticated diagnostic tools, the decisive factor underlying this progress is the parallel ongoing development of the entire range of dental instruments.

With its new and innovative products, Aesculap is playing a decisive part in this development. Aesculap's products are developed strictly by the close contact to our customers.

This applies as much to new and innovative areas such as micro-surgical instruments for implantology, periodontics and endodontics, as to standard dental instruments in general. Throughout the world, Aesculap products have proved themselves in day-to-day practice, not least from the point of view of quality and cost. With limited investment budgets available,

optimized hospital procedures and resulting efficiency gains offer a means of making better use of existing resources. Aesculap's range of services can help to achieve such improvements.

The following service components enable us to meet the ever increasing demands of the future:

- our many years of experience and competence in the field of planning and advice,
- the world renowned excellent quality of our products,
- our high level of product availability,
- individual and practice-oriented service concepts to maintain the value of your dental instruments,
- continuous new opportunities for ongoing and further training in the form of the courses offered by the Aesculap Academy.

With ERGOPLANT, the ergonomically designed instrument for dental implantology, Aesculap is setting new standards:

For the first time a complete and unitary product-line, which covers all indications, is available to the implantology focused dentist. Where ever necessary, the working tips had been optimized for each indication. Due to the ergonomic handle design the user is supported by easiness and grip. Finally the surgeon can concentrate on the patient respectively the surgery. The instrument handling is done intuitively.

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
In close collaboration with Dr. Dr. Andres Stricker of Freiburg University, an instrument set especially designed for dental implantology has been developed for optimal ergonomics concerning the working tips and handle for any specific indication.

ERGOPLANT thus enables working efficiently with minimum physical strain under practice conditions.

All ERGOPLANT instruments are suitable for steam sterilization 134°C (Autoclav).

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Implantology Basic Kit

The standard instrumentation for all dento-alveolar operations.

DX705

Aesculap Basic Kit for dental implantology, consisting of:

BB515	1 PAK = Package of 100 scalpel blades, sterile, Fig. 15
BB063R	1 Scalpel handle, No. 3
DX094	1 ERGOPLANT Mouth mirror holder with PEEK handle
DA036R	1 PAK = Package of 12 Rhodium-Mouth mirror, 22 mm diam.
DX051R	1 ERGOPLANT Tissue forceps
DX052R	1 ERGOPLANT Dissecting forceps
DX201R	1 ERGOPLANT Raspatory, BUSER
DX200R	1 ERGOPLANT Raspatory, PRICHARD
DX220R	1 ERGOPLANT Bone curette, LUCAS
BM065R	1 DUROGRIP Needle holder, HEGAR-MAYO
DO219R	1 Gum scissors, LOCKLIN
DX500R	1 ERGOPLANT Bone rongeur, LUER-FRIEDMANN
JF244R	1 Tray (274 x 172 x 30 mm) suitable for dental container JN092
TE836	1 Packing stencil (only available for DX705)

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Implantology Advanced Kit
Special instrumentation designed for the particular requirements of the bony implant site.

With this instrumentation it is possible to remove inflamed soft tissue residue completely while achieving maximum preservation of alveolar bone.

This reduces the need for augmentative measures.

ERGOPLANT – Advanced Bone Preparation

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
DX706

Ideal for detailed preparation of the bony implant site and for thorough removal of soft tissue residue. Also suitable for minor augmentative procedures.

Aesculap Advanced Kit for dental implantology, consisting of:

DX230R	1 ERGOPLANT Excavator
DX231R	1 ERGOPLANT Excavator
DX241R	1 ERGOPLANT Universal curette, M23A
DX300R	1 ERGOPLANT Furcation probe, NABERS
DX400R	1 ERGOPLANT Periotome
DX253R	1 ERGOPLANT Bone file, MILLER
DX252R	1 ERGOPLANT Bone file, MILLER
DX555R	1 ERGOPLANT Applicator for bone material, STRICKER
DX600R	1 ERGOPLANT Membrane placing instrument, STRICKER
DX501R	1 ERGOPLANT Bone cutting forceps
JF244R	1 Tray (274 x 172 x 30 mm) suitable for dental container JN092
TE839	1 Packing stencil (only available for DX706)

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
Filigree-aesthetic soft tissue management in implantological operations as the precondition for a good implantation.

These instruments can be used as successfully for larger or coarser mucocutaneous incisions as for the sensitive, easily torn mucous membrane in the area of the neighbouring teeth and their papillae. In soft tissue augmentations, such as free mucous membrane or connective tissue transplants, both harvest and transplantation site can be prepared simultaneously in exemplary fashion.

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DX707

Aesculap Soft Tissue Kit for dental implantology, consisting of:

BB045T	1 Titanium scalpel handle
BB367R	1 PAK = Package of 10 micro scalpel blades, sterile, Fig. 367
BB063R	1 Scalpel handle, No. 3
BB515	1 PAK = Package of 100 scalpel blades, sterile, Fig. 15
DX050R	1 ERGOPLANT Tissue forceps
BM003R	1 DUROGRIP Needle holder
DX450R	1 ERGOPLANT Gingivectomy knife, KIRKLAND
DX203R	1 ERGOPLANT Raspatory
DX202R	1 ERGOPLANT Raspatory, STRICKER
DX405R	1 ERGOPLANT Atraumatic papilla elevator, PHW
DX401R	1 ERGOPLANT Spatula, sharp, STRICKER
DX402R	1 ERGOPLANT Spatula, HEIDEMANN
BC061R	1 Delicate scissors, curved
DX055R	1 ERGOPLANT Soft Tissue Forceps
JF244R	1 Tray (274 x 172 x 30 mm) suitable for dental container JN092
TE837	1 Packing stencil (only available for DX707)

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
In case of poor bone thickness up to 3mm, it's possible to split the jaw-bone (Bone Split) and separate it by means of special chisels, the so called Osteotomes. Finally a gap occurs between buccal and lingual respectively buccal and palatinal bone-lamella, where we can insert an implant.

Another indication could be the "chiseling out" of bone-transplants. In any case the user is supported by the ERGOPLANT handle design. Together with Dr. Stricker, we developed a special "Step-Osteotome". Due to this innovative instrument, cracks on the bone-lamella, as often occurs particularly on thin jaw-bones, can be avoided now.

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DX708

Aesculap Bone Split Kit for dental implantology, consisting of:

DX510R	1 ERGOPLANT Osteotome, straight, 3 mm
DX511R	1 ERGOPLANT Osteotome, straight, 4 mm
DX512R	1 ERGOPLANT Osteotome, straight, 6 mm
DX513R	1 ERGOPLANT Osteotome, straight, 8 mm
DX517R	1 ERGOPLANT Osteotome, curved, 4 mm
DX518R	1 ERGOPLANT Osteotome, curved, 6 mm
DX522R	1 ERGOPLANT Step osteotome, 6 mm
JF283R	1 Tray (274 x 172 x 30 mm) suitable for dental container JN092
TE885	1 Packing stencil (only available for DX708)

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
Under circumstances the achievement of primary stability in implantations could be problematically. For instance if the upper jaw's bone structure is soft and "comb-like". With the Bone Condenser Kit the bone can be condensed (Bone Condensing) step-by-step until we finally achieve enough primary stability. The aligned conical working tips are convex on the tip. The rather spongy bone is almost "transformed" into compact cortical bone, which now enables a secure and successful implant insertion.

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DX709

Aesculap Bone Condenser Kit for dental implantology, consisting of:

DX530R	1 ERGOPLANT Bone condenser, straight, 1.5 - 2.2 mm
DX531R	1 ERGOPLANT Bone condenser, straight, 2.2 - 2.8 mm
DX532R	1 ERGOPLANT Bone condenser, straight, 2.8 - 3.5 mm
DX533R	1 ERGOPLANT Bone condenser, straight, 3.5 - 4.2 mm
DX534R	1 ERGOPLANT Bone condenser, straight, 4.2 - 4.8 mm
DX535R	1 ERGOPLANT Bone condenser, bayonet, 1.5 - 2.2 mm
DX536R	1 ERGOPLANT Bone condenser, bayonet, 2.2 - 2.8 mm
DX537R	1 ERGOPLANT Bone condenser, bayonet, 2.8 - 3.5 mm
DX538R	1 ERGOPLANT Bone condenser, bayonet, 3.5 - 4.2 mm
DX539R	1 ERGOPLANT Bone condenser, bayonet, 4.2 - 4.8 mm
JF283R	1 Tray (274 x 172 x 30 mm) suitable for dental container JN092
TE886	1 Packing stencil (only available for DX709)

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
In case of upper-jaw augmentation, it is possible to pass the external sinus elevation with its lateral window technique, if still enough vertical bone is available. Then an internal Sinus Elevation by means of the Sinutome Kit can be proceeded. With different sizes of Sinutomes the maxillary sinus floor is elevated through the oral cavity towards cranial. The concave working tip is semi-sharp on the edge, which enables the taking along of additional bone for further elevation. After placing the augmentation-material, the implant can be inserted in the same session.

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DX710

Aesculap Sinutome Kit for dental implantology, consisting of:

DX570R	1 ERGOPLANT Sinutome, straight, Ø 2.2 mm
DX571R	1 ERGOPLANT Sinutome, straight, Ø 2.8 mm
DX572R	1 ERGOPLANT Sinutome, straight, Ø 3.5 mm
DX573R	1 ERGOPLANT Sinutome, straight, Ø 4.2 mm
DX574R	1 ERGOPLANT Sinutome, straight, Ø 4.8 mm
DX575R	1 ERGOPLANT Sinutome, bayonet, Ø 2.2 mm
DX576R	1 ERGOPLANT Sinutome, bayonet, Ø 2.8 mm
DX577R	1 ERGOPLANT Sinutome, bayonet, Ø 3.5 mm
DX578R	1 ERGOPLANT Sinutome, bayonet, Ø 4.2 mm
DX579R	1 ERGOPLANT Sinutome, bayonet, Ø 4.8 mm
JF283R	1 Tray (274 x 172 x 30 mm) suitable for dental container JN092
TE887	1 Packing stencil (only available for DX710)

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
Autogenous bone in augmentation has been established as the "golden standard". The region where to gain the bone (Bone Harvesting) and the required consistence for augmentation can vary, depending on the case. Here the ERGPLANT Micro Bone Mill shows its flexibility whereat the maximum yield of bone-grist is always ensured.

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DX800

Aesculap Micro Bone Mill Kit,
consisting of:

- | | |
|--------|----------------------------------------------------------------------------------------------------------------|
| DX801R | 1 ERGOPLANT Micro bone-mill
(incl. cutting disk, fine) |
| DX803R | 1 ERGOPLANT Cutting disk, coarse |
| DX810R | 1 ERGOPLANT Bone mill-cleaner, fine |
| DX811R | 1 ERGOPLANT Bone mill-cleaner, coarse |
| JF284R | 1 Siebtray (274 x 172 x 60 mm)
with storage aid for micro bone-mill,
suitable for dental container JN095 |
| TE888 | 1 Packing stencil (only available for DX800) |

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

Delivered with fine cutting disk
DX802R

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ERGOPLANT Micro bone-mill

- small and handy
- interchangeable cutting disks
- cutting disks available for finer or coarser bone particles
- easy to clean due special cleaning instruments
- easy mounting and demounting
- maximum bone yield requires less bone harvesting

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Recommended accessories for the micro bone-mill:
Cleaning instruments DX810R and DX811R as shown on page 20.

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Cleaning spatula for radial bone material collecting by cleaning the cutting disk's edge (DX811R the complete surface) at the same time. Valuable bone material can be used for augmentation.

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Cleaning fork for punching out valuable bone material of all cutting disk holes at once in a row.

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ERGOPLANT
Bone-mill-cleaner fine
for DX802R
195 mm, 7³/₄"

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ERGOPLANT
Bone-mill-cleaner coarse
for DX803R
195 mm, 7³/₄"

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
Enables secure carrying and placing of enough bone material. The other working tip serves for pushing and plugging bone material into the augmentation area.

**STRICKER
DX555R**
ERGOPLANT
Aufbringinstrument
für Knochenmaterial
ERGOPLANT
Applicator for
bone material
195 mm, 7 ³/₄"

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
The small surface of the pointed tip avoids the sticking of the humid membrane. The small spatula on the other end, helps to take out bone material of, for instance, the bone files DX252R-DX253R, as shown on page 23.

**STRICKER
DX600R**
ERGOPLANT
Membran
Platzierungsinstrument
ERGOPLANT
Membrane placing
instrument
190 mm, 7 ¹/₂"

LUER-FRIEDMANN
DX500R
 ERGOPLANT
 Hohlmeißelzange
 ERGOPLANT
 Bone rongeur
 145 mm, 5¾"

DX501R
 ERGOPLANT
 Knochensplitterzange, gebogen
 ERGOPLANT
 Bone cutting forceps, curved
 130 mm, 5"

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For detailed grabbing and removing
 of granulated tissue. Also suitable for
 cutting smaller bone transplants.

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

If smaller quantity of bone material is required, with the bone files DX252R-DX253R it can be harvested directly from the jaw-bone, without additional trauma for the patient. The quantity is often sufficient, for augmentation around the implant-neck.

ERGOPLANT Immediate Implant Retractor

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The ERGOPLANT Immediate Implant Retractor ensures a safe implantation by holding back the alveolar soft tissue.

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Safe implantation without contamination

In case of immediate implantation there is a great risk to contaminate the jaw-bone with connective tissue fibre of the soft tissue. The ERGOPLANT Immediate Implant Retractor ensures a safe implantation by holding back the alveolar soft tissue. Thereby working will be safer, easier and more comfortable.

DX110R

ERGOPLANT Immediate
Implant Retractor
ERGOPLANT Immediate
Implant Retractor

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
ERGOPLANT Mallet 135 g with interchangeable plastic disks, 200 mm, 8"

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
The Step-Osteotome was specially developed for secure bone-splitting with simultaneous stepwise widening. Due to the specially designed working tips, the maxillary crest will be first expanded towards palatinal and only afterwards towards vestibular. In the time spongy bone in the center will be condensed. After removing the Step-Osteotome, a pre-selection point for the implant bur can be clearly seen.

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

The maxillary crest is splitted transversal by the step osteotom.

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
After the bone split is kept apart by the retractor.

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The implant can be positioned optimally by using the existing bone.

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After splitting the maxillary crest the bone split is closing within a short period of time. This leads to an unwished and unnecessary bone loss when drilling afterwards. This causes the risk of an implant that is not ideally positioned and additional augmentative treatments could be necessary.

The bone split is kept apart by the retractor. Thus, the existing bone can be used optimally to position the implant.

ERGOPLANT – Bone Split Retraction System

DX505R

ERGOPLANT
Bone Split Retraktor-
Applizierinstrument

ERGOPLANT
Applicator for
Bone split retractors

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





Separating the maxillary sinus membrane carefully from the bony base.
Finally the sinus window will be pushed with the other working tip gently into the hollow maxillary cave.

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The Schneiderian membrane of the sinus can be detached with elevator D0571R by 1- 2 mm along the total cortical bone

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Elevator D0573R is particularly useful in the lower part for detaching the membrane further. The same goes for the caudal area, especially when the membrane is 4-5 mm above the base of the maxillary sinus and the maxillary sinus is built up with granular material.

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Elevator D0575R makes it possible to lift the membrane further in mesial region. One working tip is specially designed for the right maxillary sinus and the other for the left maxillary sinus.

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With elevator D0574R, the sinus membrane is further elevated in the upper and posterior part of the maxillary sinus wall.

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

With the aid of elevator D0572R, the sinus membrane is further detached over a broad area in the inner region, especially along the roof of the palate. Condenser D0580R permits positioning of the filling material.

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Reduced sinus lift instrument set.

Two double-ended sinus lift instruments consisting of D0576R (5 mm condenser/spoon elevator) and D0577R (90° dissector / raspatory).

Many years of experience with sinus lift instruments has shown that this reduced instrument set is more than sufficient for surgical procedures.

The reduced set appreciably lowers the investment the hospital or clinic needs to make.

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
This illustration shows the technique
of the percrestale sinus elevation.

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
ERGOPLANT
Tissue forceps,
with plateau
120 mm, 4 $\frac{3}{4}$ "

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ERGOPLANT
Non-traumatic
dissecting forceps
120 mm, 4 $\frac{3}{4}$ "


$\frac{1}{1}$

DX055R

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ERGOPLANT
Soft tissue forceps
140 mm, 5 $\frac{1}{2}$ "

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
ERGOPLANT
Spatula, sharp
special for soft tissue
190 mm, 7 1/2"

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


Atraumatic dissecting of the flap by
means of the papilla elevator without
any damage to the papillae.

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



DX720

Micro Soft Tissue Set acc. Steigmann,
consisting of:

- DX305R Handle for micro blades, 150 mm
- DX301R Micro raspatory, 180 mm
- DX303R Micro tissue forceps, 180 mm
- FD252R Micro ligature scissors, curved,
one blade serrated, 180 mm
- FD253R Micro dissecting forceps, 180 mm
- BM003R DUROGRIP Needle holder, 170 mm
- JG381R Tray, 274 x 172 x 30 mm

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

Micro ligature scissors, curved,
one blade serrated

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Micro dissecting forceps,
diamond dust coated

180 mm, 7"

ERGOPLANT – Soft Tissue Surgery

JG381R

Siebtray, 274 x 172 x 30 mm,
passend in Dental Container JN092

Tray, 274 x 172 x 30 mm,
suitable for dental container JN092

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
This working tip helps in case of suturing mucosa transplants, or as raspatory for larger surfaces.

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
ERGOPLANT
Raspatory for fine soft tissue
190 mm, 7 1/2"

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
ERGOPLANT
Raspatory for soft tissue
sharp/blunt
195 mm, 7 3/4"

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


$\frac{1}{1}$

DX400R
ERGOPLANT
Periotom
ERGOPLANT
Periotome
170 mm, 6 $\frac{3}{4}$ "


ERGOPLANT – Advanced Bone Preparation

SCHLUGER
DX250R
FS9/10S
ERGOPLANT
Feile
ERGOPLANT
File
185 mm, 7 1/4"

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
To ensure a successful long-term implantation, a consequent post-operative treatment is very important. To eliminate plaque on the implant – abutment – interface without damaging the implant titanium curettes are the optimal choice.

ERGOPLANT – Recall Titanium Currettes

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These instruments were developed to hold back lips, cheeks, tongue and lingual mucous membrane. Due to the ergonomic handle design, a more relaxed and fatigue-free working is possible, even in long operations.

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

Sales unit: PAK = Package of 2 pieces
suitable for steam sterilization 134°C
(Autoclave)

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The mouth wedges are particu-
larly advantageous for longer
implantations and also provide
support for the patient.

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ERGOPLANT
Silicone mouth wedge
Size: S

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ERGOPLANT
Silicone mouth wedge with metal core
Size: M

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ERGOPLANT
Silicone mouth wedge with metal core
Size: L

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

The start of a new extraction culture.


With a new method to extract teeth in the shortest possible time without damaging the alveoli. The innovative AESCULAP exCenter extraction system offers safe and modern tooth extraction with many advantages for both user and patient.

Advantages over conventional methods:

- Almost no trauma for the patient
- Sensitive extraction with no force required
- No lever pressure on tooth
- Practically eliminates tooth crushing
- Alveoli remain to a large extent undamaged
- Preserves the buccal lamella
- Also ideal for immediate implantation
- Prevents inflammation and abscess formation
- Only one extraction forceps for all teeth
- Conventional tooth extraction forceps not required
- Ergonomic design
- Clearly organised instrumentation
- Safe preparation and storage in AESCULAP trays

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exCenter II extraction forceps without working tip,
universal for upper- and lower jaw,
function only with working tips DH822R - DH825R
(to be ordered separately)

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exCenter II working tip 1
for upper incisors, upper and
lower Canines, Premolars

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
exCenter II working tip 2
for lower incisors, upper and
lower roots

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exCenter II working tip 3
for upper and lower Molars,
on intact crowns

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exCenter II working tip 4
"Bifurcation"
for upper and lower Molars,
on partly damaged crowns

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
DH820

Aesculap exCenter II Extraction System,
consisting of:

- | | |
|--------|-------------------------------------------------------------------|
| DH821R | 1 Eccentric II tooth forceps without tip |
| D0801R | 1 Desmotome, straight, filigree |
| D0802R | 1 Desmotome, left, filigree |
| D0803R | 1 Desmotome, right, filigree |
| JF267R | 1 Tray (274 x 172 x 30 mm)
suitable for dental container JN092 |
| TE884 | 1 Packing stencil (only available for DH820) |

Sterilgutkreislauf – Instrumentenaufbereitung Sterile Goods Cycle – Instrument Reprocessing

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
The sterile goods cycle in the dental practice with the Aesculap Sterile Container System

All hygienic measures in the dental practice are primarily for keeping healthy the patients, the dentist and the assistants.

To ensure the best possible safety it is particularly important to observe the principles of hygiene when reprocessing instruments.

To keep practice work safe, efficient and economical, clear procedures must be specified for practice hygiene.

The Aesculap Sterile Container System has proven itself as an excellent aid for organizing the instrument reprocessing procedure.

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Art.Nr. Cat.No.	Seite Page	Art.Nr. Cat.No.	Seite Page	Art.Nr. Cat.No.	Seite Page	Art.Nr. Cat.No.	Seite Page
BB045T	9, 58	DO804R	64	DX402R	9, 42	DX579R	5, 39
BB063R	5, 9, 59	DO805R	64	DX405R	9, 42	DX600R	7, 21
BB363R	58			DX450R	9, 40	DX705	5
BB364R	58	DX050R	9, 41	DX500R	5, 22	DX706	7
BB367R	9, 58	DX051R	5, 41	DX501R	7, 22	DX707	9
BB369R	58	DX052R	5, 41	DX505R	28, 29	DX709	13
BB511	59	DX055R	9, 41	DX506T	29	DX710	15
BB512	59	DX094	5, 40	DX507T	29	DX720	44
BB515	5, 9, 59	DX100R	54	DX510R	11, 26	DX800	17
BB542	59	DX101R	55	DX511R	11, 26	DX801R	17, 18, 19
		DX102R	55	DX512R	11, 26	DX802R	18, 19
BC061R	9, 60	DX105R	56	DX513R	11, 27	DX803R	17, 18, 19
		DX106R	56	DX517R	11, 27	DX807	11
BH111R	60	DX110R	22	DX518R	11, 27	DX810R	17, 20
		DX113	57	DX522R	11, 25	DX811R	17, 20
BM003R	9, 44, 47, 61	DX114	57	DX530R	13, 30, 31		
BM065R	5, 61	DX115	57	DX531R	13, 30, 31	FD252R	44, 46
		DX120R	57	DX532R	13, 30, 31	FD253R	44, 46
DA036R	5, 40	DX200R	5, 43	DX533R	13, 30, 31		
DA037R	40	DX201R	5, 43	DX534R	13, 30, 31	GD580R	59
		DX202R	9, 48	DX535R	13, 32	GD584R	59
DH820	65	DX203R	9, 48	DX536R	13, 32	GD585R	59
DH821R	63, 65	DX220R	5, 49	DX537R	13, 32		
DH822R	63	DX230R	7, 49	DX538R	13, 32	JF244R	5, 7, 9
DH823R	63	DX231R	7, 49	DX539R	13, 32	JF267R	65
DH824R	63	DX232R	49	DX545R	25	JF283R	11, 13, 15
DH825R	63	DX240R	50	DX550R	23	JF284R	17
		DX241R	7, 50	DX555R	7, 21		
DO219R	5, 60	DX250R	51	DX560R	33	JG381R	44, 47
DO571R	34	DX252R	7, 23, 51	DX561R	33		
DO572R	36	DX253R	7, 23, 51	DX570R	15, 38	TE836	5
DO573R	34	DX300R	7, 51	DX571R	15, 38	TE837	9
DO574R	35	DX301R	43, 44, 45	DX572R	15, 38	TE839	7
DO575R	35	DX303R	44, 45	DX573R	15, 38	TE884	65
DO576R	36	DX305R	44, 45, 58	DX574R	15, 38	TE885	11
DO577R	37	DX360T	53	DX575R	15, 39	TE886	13
DO801R	64, 65	DX361T	53	DX576R	5, 39	TE887	15
DO802R	64, 65	DX400R	7, 50	DX577R	5, 39	TE888	17
DO803R	64, 65	DX401R	9, 42	DX578R	5, 39		

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Aesculap AG & Co. KG
B. Braun Austria GmbH