www.maslanka.de

H.+ H. Maslanka  
Chirurgische Instrumente GmbH  
Stockacher Str. 172  
D-78532 Tuttlingen

CE 0483

☎ 0049-(0)7461-9607-0  
📠 0049-(0)7461-4732  
✉ info@maslanka.de

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**Seite 2 - 8**

**Flex-Pusher II (Einweg)**

**EN Instruction for Use**

**Pages 9 – 15**

**Flex-Pusher II (disposable)**

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**Pagina 16 – 22**

**Flex-Pusher II (monouso)**

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**Page 23 – 29**

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**Página 30 – 38**

**Flex-Pusher II (desechable)**

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**Páginas 37 – 45**

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**Strona 44 - 50**

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### Flex-Pusher II (disposable)

---

Thank you for choosing one of the high-quality products from our company. For the safest possible use of the product, we have summarised all important information in this guide. We are grateful for additional advice and suggestions.

#### Warning – Important information

Among other things, complications may occur with the application of the Flex-Pusher.

- Cardiac arrhythmia
- Formation of blood clots, vascular injury and occlusion due to existing vascular calcification and heart disease
- Unintentional injuries to the heart
- Injuries to the blood vessels
- Bleeding or secondary bleeding in the heart or pericardium due to injury of the heart
- Infection via the point of introduction
- Pulmonary embolism
- Arterial embolism
- Cardiac tamponades
- Perforation
- Apnea
- Arrhythmias
- Attak (TIA) or stroke
- Vomiting
- Fever
- Heart block
- Cardiac arrest
- Hypotension
- Infection including endocarditis
- Headache / migraine
- Myocardial infraction
- Need for surgery
- Perforation of the vessel or myocardium
- Anesthesia reaction
- Subaortic stenosis
- Death
- AV-fistula
- Air embolism
- Pericardial effusion
- Pericarditis
- Post-pericardiotomy Syndrome
- Pseudoaneurysm
- Pulmonary edema
- Stroke

Appropriate preparations for any ensuing complications must be taken before starting the application!

#### Warning – Important information

- The Flex-Pusher may only be used in connection with compatible heart implants from Occlutech!
- Please observe the colour code; the coil colour of the Flex-Pusher must agree with the label colour of the occluder!
- The Flex-Pushers or occluders may not be interchanged!
- Always use only one unit!
- Please ensure that you follow the instructions for use for the Occlutech closure system and introduction system!

#### Warning – Important information

Should the occluder loosen prematurely for some unknown reason, it is possible to retrieve the occluder with alligator forceps or a sling and a larger sheath.

An emergency kit for the retrieval of the occluder must be available in the room!

Due to the possibility of serious complications occurring, a cardiac machine with appropriate staff for an emergency operation must be immediately available on site!

## 1. Intended use

The Flex-Pusher is an accessory compatible heart implants from Occlutech and is applied with special cardiac catheters in the human heart. Alongside it, the jointly applicable instructions for use for Occlutech implants must also be taken into consideration.

The Flex-Pusher positions the Occlutech heart implants on the atrium septum, ventricle septum or on the left atrial appendage via a minimally invasive catheter introduction system.

- These instruments are not for other applications!
- Damaged instruments should not be used!

These instructions for use should help the user to use the Flex-Pusher for the intended purpose, to convey how they are to be handled and through this, ensure greater safety for the patient.

### 1.1 Indication

Flex-Pusher is used to position the Occlutech heart implants, which are implants used for closure of defects in septum of the heart and the left atrial appendage.

### 1.2 Contraindication

The contraindications are Occlutech heart implants related and listed in the corresponding instruction for use.

## 2. Designs

**H. + H. Maslanka, Chirurgische Instrumente GmbH** offers the Flex-Pusher in various designs which shall be used with Occlutech's colour-corresponding products. For information regarding the dimensions which are available please ask your supplier, Occlutech, or ask us directly.

Only use the Flex-Pusher and occluder with the introducer sheath designed for that purpose (see table below).

Article number Flex-Pusher	Article number Occluder	Instrument diameter	Introducer sheath internal diameter
51FP060	29 ASD 05, 29 ASD 04	1,90mm	7F
	71 VSD 04	1,90mm	6F
51FP100	29 ASD 10, 29 ASD 09, 29 ASD 07, 29 ASD 06, 19 PFO 18D, 16 UNI 17, 71 VSD 08	2,15mm	7F
	71 VSD 06	2,15mm	6F
	71 VSD 10	2,15mm	8F
51FP120	29 ASD 18, 29 ASD 15, 29 ASD 12, 29 ASD 13, 29 ASD 16, 18 PFO 25S, 19 PFO 30D, 19 PFO 25D, 16 UNI 28, 16 UNI 24, 71 VSD 16	2,65mm	9F
	71 VSD 14, 71 VSD 12	2,65mm	8F
	71 VSD 18	2,65mm	10F
	71 VSD 20	2,65mm	11F
	30 LAA 18, 30 LAA 21	2,65mm	12F
	30 LAA 24	2,65mm	14F
51FP150	29 ASD 24, 29 ASD 21	3,26mm	11F
	30 LAA 27	3,26mm	14F
51FP160	29 ASD 40, 29 ASD 39, 29 ASD 36, 29 ASD 33, 29 ASD 30, 29 ASD 27, 16 UNI 40	3,26mm	12F
	19 PFO 35D, 16 UNI 33	3,26mm	11F
	30 LAA 30, 30 LAA 39	3,26mm	14F

## 3. Before the application

The user should carefully study these instructions for use as well as Occlutech heart implants and catheter before using the instrument. He needs to have fully understood how the instrument is operated and handled.

After removing the instrument from the packaging, it must be checked for defects and damage; this includes the expiration date of the product. Should damage or defects be identified whilst checking, the instrument should not be used on patients and should be promptly returned to **H. + H. Maslanka, Chirurgische Instrumente GmbH** for evaluation.

The wearing of protective clothing (goggles, gloves, surgical clothes, etc.) is essential.

**Flex-Pushers are supplied sterile for single use. Prior to the application the instrument must be checked for damage. Only intact and sterile instruments should be used! The instrument should only be used once. Please compare the size of the Flex-Pusher with the size of the catheter!**

## 4. Handling

### 4.1 Pre-bending of the Flex-Pusher (only applies to shapeable design)


Flex-Pushers with a bendable spring can be pre-bent at the distal end before use. (See article description on packaging.) There should be no additional bends.

**The Flex-Pusher must always be bent from the distal end to the handle! Do not use instruments that are damaged during the bending process! Bent Flex-Pushers may not be bent back.**

### 4.2 Preparation


1. Insert the Flex-Pusher with the distal end first through the (fitting) introducer.

2. Then couple the supplied occluder to the Flex-Pusher.

For this, the pick up (hook) is moved out by actuating the coil on the handle and the adapter ball of the occluder is inserted into the pick up.

3. By actuating or letting the coil go, you re-retract the pick up. Turn the locking nut as far as the coil in order to prevent the unintentional de-coupling of the occluder from the Flex-Pusher

4. Please ensure that the occluder is firmly and securely coupled to the Flex-Pusher. Turn the locking nut as far as it will go (preferably in a straight position).

**Carry out the following steps steps in a sterile, physiological NaCl solution**

5. Now fully pull in the occluder into the introducer. This is done by simultaneously holding the introducer tight and pulling back the Flex-Pusher.


**This must take place in a physiological NaCl solution in order to avoid introducing air bubbles into the introducer.**

6. Then, through the rinse input, re-rinse the introducer with the occluder located within it using the NaCl solution and check it for any air bubbles.

### 4.3 Introducing the catheter (PFO, ASD and UNI Occluder)

1. Please make sure that there is no air in the long delivery sheath (catheter) and the introducer, as at this point, there is a risk of air embolism. The introducer is now inserted into the catheter. To avoid kinking the long delivery sheath, it is necessary to close the occluder flush with the end of the long delivery sheath.
2. Connect the introducer by using the screw connector on the delivery sheath.
3. Now, together with the occluder, slowly slide the Flex-Pusher from the introducer into the long delivery sheath until its tip (The tip of the catheter is located in the left atrium).
4. Then, by carefully pushing the Flex-Pusher further, the left atrial disc of the occluder is opened.

5. Hold the occluder in this position and pull back carefully the entire introduction system (catheter and Flex-Pusher) against the septum until an elastic resistance with pulse-synchronous movements is established and in TEE (transesophageal echocardiography), the left atrial disc abut.
6. Constantly pulling on the catheter, the right atrial disc is slowly opened in the right atrium until it is fully unfolded.

7. Through the application of a contrast agent in the right atrium establish whether the occluder has been correctly positioned in order to rule out any residual shunt or situations where the occluder does not optimally abut.
8. If it is established that the occluder is not correctly positioned or is too large and/or small, the occluder can be refed by reversing points 4.3.6 to 4.3.3.  
**At this stage, particular care is required!**
9. After positioning of the occluder has taken place, you can release the occluder from the Flex-Pusher by turning back the locking nut and slowly pushing the coil forward. Please do not press the Flex-Pusher against the occluder in the process; instead, pull it back slightly.
10. The occluder de-couples itself from the Flex-Pusher.

**Due to the anatomy of the heart, circumstances may arise where the geometric arrangement of the occluder and the Flex-Pusher is unfavourable and as such, the decoupling is made more difficult. In order to solve the problem, slide the Flex-Pusher forward slightly and attempt to release by turning the fixture between the Flex-Pusher and the occluder. Make sure that the hook of the Flex-Pusher does not penetrate into the cardiac tissue/occluder tissue and only carry out this procedure by x-ray.**



11. Remove the entire introduction system (catheter and Flex-Pusher) from the patient. Make sure that the Flex-Pusher is located within the lock. Otherwise, there is an increased risk of injury!


**Should the occluder loosen prematurely for some unknown reason, it is possible to retrieve the occluder with alligator forceps or a sling and a larger sheath!!!**

#### 4.4 Introducing the catheter (LAA Occluder)

1. Please make sure that there is no air in the long delivery sheath (catheter) and the introducer, as at this point, there is a risk of air embolism. The introducer is now inserted into the catheter. To avoid kinking the long delivery sheath, it is necessary to close the occluder flush with the end of the long delivery sheath.
2. Connect the introducer by using the screw connector on the delivery sheath.
3. Now, together with the occluder, slowly slide the Flex-Pusher from the introducer into the long delivery sheath until its tip (The tip of the catheter is located in the LAA).
4. After confirmation of correct position by TEE and fluoroscopy by carefully pushing the Flex-Pusher further, the loops are slowly pushed out from the delivery sheath

5. When unfolded a visual inspection should be performed to secure that the loops are folded in the correct position. If they are not, the occluder should be retracted in to the sheath and new attempt should be made.
6. Once correct positioning of the loops has been confirmed the body of the device is pushed out. Correct position should be controlled and a tug test shall be performed to show that the occluder is securely anchored in the LAA.

7. If the device passes the tug test, it shall be released
8. The occluder is released from the Flex-Pusher by turning back the locking nut and slowly pushing the coil forward. Please do not press the Flex-Pusher against the occluder in the process; instead, pull it back slightly
9. The occluder de-couples itself from the Flex-Pusher.

*The locking nut when open*

**Due to the anatomy of the heart, circumstances may arise where the geometric arrangement of the occluder and the Flex-Pusher is unfavourable and as such, the decoupling is made more difficult. In order to solve the problem, slide the Flex-Pusher forward slightly and attempt to release by turning the fixture between the Flex-Pusher and the occluder. Make sure that the hook of the Flex-Pusher does not penetrate into the cardiac tissue/occluder tissue and only carry out this procedure by x-ray.**


10. Remove the entire introduction system (catheter and Flex-Pusher) from the patient. Make sure that the Flex-Pusher is located within the lock. Otherwise, there is an increased risk of injury!

**Should the occluder loosen prematurely for some unknown reason, it is possible to retrieve the occluder with alligator forceps or a sling and a larger sheath!!!**

#### 4.5 Introducing the catheter (mVSD Occluder)


1. Please make sure that there is no air in the long delivery sheath (catheter) and the introducer, as at this point, there is a risk of air embolism. The introducer is now inserted into the catheter. To avoid kinking the long delivery sheath, it is necessary to close the occluder flush with the end of the long delivery sheath.
2. Connect the introducer by using the screw connector on the delivery sheath.
3. **Use TEE / TTE and angiography as a guide during each step of device deployment.**
4. Now, together with the occluder, slowly slide the Flex-Pusher from the introducer into the long delivery sheath until its tip (The tip of the catheter is located in the left ventricle).
5. By carefully pushing the Flex-Pusher further, the distal disc of the occluder is opened

6. Pull the entire assembly (catheter and Flex-Pusher) into the VSD
7. Pull the catheter further back to deploy the occluder in the defect
8. After conforming the position and under constant train on the catheter the proximal disc in the right ventricle open slowly until it unfolds completely

9. Through the application of a contrast agent in the left ventricle the correct positioning and proper closure should be ascertained. The correct positioning of the occluder has been confirmed through a tug test to show that the device is securely anchored in the septum.
10. If it is established that the occluder is not correctly positioned or is too large and/or small, the occluder can be re-fed by reversing points 4.5.6 to 4.5.3.  
**At this stage, particular care is required!**
11. The occluder is released from the Flex-Pusher by turning back the locking nut and slowly pushing the coil forward. Please do not press the Flex-Pusher against the occluder in the process; instead, pull it back slightly
12. The occluder de-couples itself from the Flex-Pusher.

*The locking nut when open*

#### 4.6 Disposal

After use, the instrument must be disposed of in accordance with statutory and internal clinical provisions

#### 5. Cleaning/Sterilisation

Flex-Pushers for single use supplied in a sterilised state using ETO gas. Repeat processing and sterilisation after the application is no longer possible!

The Flex-Pusher may only be used once.

Re-sterilisation after the shelf life should not be carried out, i.e. the product should be disposed of according to clinical provisions.

**The raw materials/materials from which the Flex-Pusher is made of can be impaired through the repeat processing and re-sterilisation of their properties.  
If the product is nevertheless reprocessed and re-sterilised, responsibility lies with the user.**

6. Notes concerning storage

For the storage duration and the shelf life, please refer to the article label.

**Do not place objects on the instruments and/or their packaging!**  
**Do not store in the vicinity of chemicals, disinfectants or radiation! Do not expose the instrument to direct or indirect sunlight or other ultraviolet radiation!**

No complaints can be made with respect to improperly-stored instruments.

7. Users

The users of **H. + H. Maslanka, Chirurgische Instrumente GmbH's** Flex-Pusher must be specialised in the field of cardiology and have sufficient experience with regard to cardiac catheterization. It is important to carefully read and observe the instructions for use. If anything is unclear, the product should not be applied. In this case, please contact **H. + H. Maslanka, Chirurgische Instrumente GmbH** or **Occlutech**.

8. Shipments relating to repairs and complaints



In the interest of the health of our employees, only disinfected and/or sterilized instruments are accepted as returns. These should be labelled as such as otherwise no further processing can take place.

A return form can be found online at:  
[www.maslanka.de](http://www.maslanka.de)

9. Legal basis

EU law is applicable.

10. Symbols used

	Not for reuse		Manufacturer
	Sterilised using ethylene oxide		Instrument diameter
	Do not re-sterilise		Warning, risks
	Store in a place which is protected from direct sunlight		Batch name
	Store in a dry place, protect from wet		Do not use if the packaging is damaged
	Observe the instruction		Use by (example)
	Instrument length		Distributed by
	Quantity unit		

11. Service

We are always here for you during our regular office hours in the case of difficulties and questions in relation to our instruments:

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