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AMPLATZER™ Membranous VSD Occluder

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AUTHOR	TITLE	JOURNAL	TYPE OF STUDY	PATIENTS ENROLLED	SUCCESSFUL DEVICE IMPLANTATION	MAJOR COMPL.	MINOR COMPL.	CLOSURE SUCCESS IN FU
S. Leng Kang	Longer-Term Outcome of periventricular device closure of muscular ventricular septal	Catheterization and Cardiovascular Interventions, 2015	Retrospective	10	100%	1	NA	90%
J. Sathananthan	Evolution in the Management of Postinfarct Ventricular Septal Defects From Surgical to Percutaneous Approach: A Single-Center Experience	Journal of Invasive Cardiology, 2013	Prospective	25	100%	NA	NA	56%
G. Butera	Over-the-Wire-Technique Device Implantation	Catheterization and Cardiovascular Interventions, 2012	Case report	5	100%	NA	NA	100%
P. Chungsomprasing	The Results of Transcatheter Closure of VSD Using Amplatzer Device and Nit Occlud Le Coil	Catheterization and Cardiovascular Interventions, 2011	Retrospective	116	100%	NA	NA	89.5%
I. Michel-Behnke	Device Closure of Ventricular Septal Defects by Hybrid Procedures: A Multicenter Retrospective Study	Catheterization and Cardiovascular Interventions, 2011	Retrospective	26	88.5%	3	NA	58%
J. S Dua	Transcatheter Closure of Postsurgical Residual Ventricular Septal Defects: Early and Mid-Term Results	Catheterization and Cardiovascular Interventions, 2009	Prospective	170	100%	NA	NA	100%
D. Kenny	Evolution of Transcatheter Closure of Perimembranous Ventricular Septal Defects in a Single Centre	Catheterization and Cardiovascular Interventions, 2009	Prospective	25	92%	2	NA	93%
H. Thiele	Immediate primary transcatheter closure of postinfarction ventricular septal defects	European Heart Journal, 2009	Prospective	29	86%	41%	NA	
D. S. Crossland	Initial Results of Primary Device Closure of Large Muscular Ventricular Septal Defects in Early Infancy Using Periventricular Access	Catheterization and Cardiovascular Interventions, 2008	Observational	8	100%	NA	NA	87.5%

Q. Zheng	A Comparative Study: Early Results and Complications of Percutaneous and Surgical Closure of Ventricular Septal Defect	Cardiology, 2008	Observational	2178	99.8%	1.9%	0.6%	99%
G. Butera	Transcatheter Closure of Perimembranous Ventricular Septal Defects - Early and Long Term Results	Journal of the American College of Cardiology, 2007	Prospective	104	96.2%	11.5%	NA	99%
M. Carminati	Transcatheter closure of congenital ventricular septal defects: results of the European Registry	European Heart Journal, 2007	Prospective	430	95%	6.5%	6.2%	95.3%
K. A. Diab	Device Closure of Muscular Ventricular Septal Defects in Infants Less Than One Year of Age Using the Amplatzer Devices: Feasibility and Outcome	Catheterization and Cardiovascular Interventions, 2007	Prospective	20	95%	0.2%	NA	100%
M. W. Martinez	Transcatheter Closure of Ischemic and Post-Traumatic Ventricular Septal Ruptures	Catheterization and Cardiovascular Interventions, 2007	Retrospective	10	100%	NA	NA	20%
C. Xunmin	Comparison of results and complications of surgical and Amplatzer device closure of perimembranous ventricular septal defects	International Journal of Cardiology, 2007	Observational	121	97%	5.5%	19%	91.8%
Y.-C. Fu	Transcatheter Closure of Perimembranous Ventricular Septal Defects Using the New Amplatzer Membranous VSD Occluder	Journal of the American College of Cardiology, 2006	Prospective	35	91%	8.6%	NA	96%
R. Holzer	Transcatheter Closure of Perimembranous Ventricular Septal Defects Using the Amplatzer Membranous VSD Occluder: Immediate and Midterm Results of an International Registry	Catheterization and Cardiovascular Interventions, 2006	Prospective	100	93%	29%	NA	83.6%
R. J. Pinto	Transcatheter Closure of Perimembranous Ventricular Septal Defects Using Amplatzer Asymmetric Ventricular Septal Defect Occluder: Preliminary Experience With 18-Month Follow up	Catheterization and Cardiovascular Interventions, 2006	Prospective	20	89.4%	NA	NA	100%
A.-A. Jameel	Retrograde Approach for Device Closure of Muscular Ventricular Septal Defects in Children and Adolescents, Using the Amplatzer Muscular Ventricular Septal Defect Occluder	Pediatric Cardiology, 2006	Case report	7	NA	NA	NA	71.4%

M. A. Walsh	Percutaneous Closure of Postoperative Ventricular Septal Defects With the Amplatzer Device	Catheterization and Cardiovascular Interventions, 2006	Case report	9	100%	0.11%	NA	55.5%
H. Akintuerk	Arterial Switch Operation for Transposition of the Great Arteries Combined With Closure of Paramembranous Ventricular Septal Defect Using the Amplatzer Membranous Ventricular Septal Occluder	American Journal of Cardiology, 2006	Case report	3	100%	NA	NA	NA
E. A. Bacha	Multicenter Experience with Periventricular Device Closure of Muscular Ventricular Septal Defects	Pediatric Cardiology, 2005	Observational	13	92.3%	0	NA	83.3%
M. Carminati	Transcatheter Closure of Congenital Ventricular Septal Defect with Amplatzer Septal Occluders	American Journal of Cardiology, 2005	Prospective	122	97.5%	NA	NA	96%
J. Masura	Percutaneous Closure of Perimembranous Ventricular Septal defects with the eccentric Amplatzer device: Multicenter follow-up study	Pediatric Cardiology, 2005	Prospective	186	100%	NA	NA	95.6%
B. D. Thanopoulos	Outcome of transcatheter closure of muscular ventricular septal defects with the Amplatzer ventricular septal defect occluder	Heart, 2005	Prospective	30	93%	3.3%	NA	93%
R. Holzer	Device Closure of Muscular Ventricular Septal Defects Using the Amplatzer Muscular Ventricular Septal Defect Occluder. Immediate and Mid-Term Results of a U.S. Registry	Journal of the American College of Cardiology, 2004	Prospective	75	86.7%	10.7%	NA	92.3%
R. Arora	Transcatheter Closure of Congenital Ventricular Septal Defects: Experience with Various Devices	Journal of Interventional Cardiology, 2003	Prospective	137	86.2% Rashkind Umbrella Device 97.1% Amplatzer VSO	NA	NA	99.2%
J. L. Bass	Initial Human Experience With the Amplatzer Perimembranous Ventricular Septal Occluder Device	Catheterization and Cardiovascular Interventions, 2003	Prospective	27	93%	NA	NA	92%
B. D. Thanopoulos	Catheter Closure of Congenital/ Acquired Muscular VSDs and Perimembranous VSDs Using the Amplatzer Devices	Journal of Interventional Cardiology, 2003	Prospective	35	95% muscular VSD 92.3% perimembranous VSD	7.6% perimembranous	NA	95%


B. D. Thanopoulos	Transcatheter closure of perimembranous ventricular septal defects with the Amplatzer asymmetric ventricular septal defect occluder: preliminary experience in children	Heart, 2003	Prospective	10	100%	1%	NA	100%
Z. M. Hijazi	Catheter Closure of Perimembranous Ventricular Septal Defects Using the New Amplatzer Membranous VSD Occluder: Initial Clinical Experience	Catheterization and Cardiovascular Interventions, 2002	Case report	6	100%	NA	NA	100%
M. Szkutnik	Postinfarction ventricular septal defect closure with Amplatzer occluders	European Journal of Cardio-thoracic Surgery, 2002	Case report	8	62.5%	NA	NA	60%
Z. M. Hijazi	Transcatheter Closure of Single Muscular Ventricular Septal Defects Using the Amplatzer Muscular VSD Occluder: Initial Results and Technical Considerations	Catheterization and Cardiovascular Interventions, 2000	Case report	8	100%	NA	NA	75%
B. D. Thanopoulos	Transcatheter Closure of Muscular Ventricular Septal Defects With the Amplatzer Ventricular Septal Defect Occluder: Initial Clinical Applications in Children	Journal of the American College of Cardiology, 1999	Case report	6	100%	NA	NA	95%

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Diametro sezione centrale (mm)	4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18
Diametro del disco destro (mm)	8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22
Diametro del disco sinistro (mm)	10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24
Materiale di Composizione	Nitinol (55% nichel 45% titanio); Dacron

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FABBRICANTE AGA Medical Corporation
USA Jude Medical Italia S.p.A.