LOTTO 48

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AMPLATZER™ PFO Occluder

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

AUTHOR	TITLE	JOURNAL	TYPE OF STUDY	PATIENTS ENROLLED	SUCCESSFUL DEVICE IMPLANTATION	MAJOR COMPL.	MINOR COMPL.	CLOSURE SUCCESS IN FU	REDUCTIO N STROKE
Patti et al	Meta-Analysis of Net Long-Term Benefit of Different Therapeutic Strategies in Patients With Cryptogenic Stroke and Patent Foramen Ovale	American Journal of Cardiology, 2015	Meta-Analysis	3311	NA	2.11%	NA	NA	68%
P. Jüni et al	Meta-analysis: Value of PFO Closure vs Medical Therapy Varies by Device	Società Italiana di Cardiologia Invasiva, 2014	Meta-Analysis	2963	NA	NA	NA	NA	77.1%


B. Meier	Patent Foramen Ovale and Closure Technique with the Amplatzer Occluder	Hindawi Publishing Corporation Scientifica, 2014	Review	308	NA	NA	NA	NA	80%
A. Migliore	Implantable devices for the closure of patent foramen ovale in adults: an Italian rapid health technology assessment	Expert Review of Medical Devices, 2014	Meta-analysis	1905	NA	1.62%	NA	NA	NA
A. Pandit	Amplatzer PFO Occluder Device may Prevent Recurrent Stroke in Patients with Patent Foramen Ovale and Cryptogenic Stroke: A Meta-Analysis of Randomised Trials	Heart, Lung and Circulation, 2014	Meta-analysis	2303	NA	0.43%	NA	NA	50%
C. A. Pickett	Percutaneous Closure versus Medical Therapy Alone for Cryptogenic Stroke Patients with a Patent Foramen Ovale: Meta-Analysis of Randomized Controlled Trials	Texas Heart Institute Journal, 2014	Meta-analysis	2303	NA	8.43%	NA	NA	67%
S. Stortecky	Percutaneous closure of patent foramen ovale in patients with cryptogenic embolism: a network meta-analysis	European Heart Journal, 2014	Meta-analysis, Randomized	2963	99.1%	62	97.8%	NA	77.1%
D. Capodanno	Updating the evidence on patent foramen ovale closure versus medical therapy in patients with cryptogenic stroke: a systematic review and comprehensive meta-analysis of 2,303 patients from three randomised trials and 2,231 patients from 11 observational studies	EuroIntervention, 2013	Meta-analysis	4534	NA	0.7%	NA	NA	77%

J. D. Carroll	Closure of Patent Foramen Ovale versus Medical Therapy after Cryptogenic Stroke (RESPECT Investigators)	The new england journal of medicine, 2013	Prospective, Randomized	980	99.1%	4.4%	96.1%	NA	75%
M. Hornung	Long-term results of a randomized trial comparing three different devices for percutaneous closure of a patent foramen ovale	European Heart Journal, 2013	Randomized	660	100%	12%	98.8%	NA	NA
G. Ntaios	PFO closure vs. medical therapy in cryptogenic stroke or transient ischemic attack: A systematic review and meta-analysis	International Journal of Cardiology, 2013	Meta-analysis	2303	NA	5.9%	NA	NA	53.2%
S. Agarwal	Meta-Analysis of Transcatheter Closure Versus Medical Therapy for Patent Foramen Ovale in Prevention of Recurrent Neurological Events After Presumed Paradoxical Embolism (MEDLINE search)	JACC: Cardiovascular Intervention, 2012	Meta-analysis, Observational	8185	99.1%	NA	NA	NA	NA
A. L. Bartorelli	Stroke Prophylaxis by Percutaneous Closure of Patent Foramen Ovale and Left Atrial Appendage	Journal of Interventional Cardiology, 2011	Review	NA	NA	NA	86%	NA	89%
N. Bissessor	Percutaneous patent foramen ovale closure: outcomes with the Premere and Amplatzer devices	Cardiovascular Revascularization Medicine, 2011	Observational	70	NA	NA	98.6%	NA	NA
A. A Khattab	Randomized clinical trial comparing percutaneous closure of patent foramen ovale (PFO) using the Amplatzer PFO Occluder with medical treatment in patients with cryptogenic embolism (PC-Trial)	Trials Journal, 2011	Prospective, Randomized	414	NA	5-8%	95%	NA	67%
G. D. Kitsios	Patent Foramen Ovale Closure and Medical Treatments for Secondary Stroke Prevention	Journal of the American Heart Association, 2011	Meta-analysis	8916	NA	NA	NA	NA	81%
C. Zhang	Transcatheter Closure of Patent Foramen Ovale in Chinese Patients With Paradoxical Embolism	Journal of the American Heart Association, 2011	Prospective	192	NA	4.16%	100%	NA	NA
D. Fischer	Transcatheter Closure of Patent Foramen Ovale (PFO) in Patients with Paradoxical Embolism: Procedural and Follow-Up Results after Implantation of the Amplatzer Occluder Device	Journal of Interventional Cardiology, 2010	Prospective	114	100%	9.7%	92%	NA	NA
A. M. Saguner	Figulla PFO Occluder versus Amplatzer PFO Occluder for Percutaneous Closure of Patent Foramen Ovale	Catheterization and Cardiovascular Intervention, 2010	Observational	40	100% AMP 95% FPO	15% FPO	100% AMP 61% FPO	NA	NA
M. Becker	Impact of Occluder Device Type on Success of Percutaneous Closure of Atrial Septal Defects-A Medium-Term Follow-up Study	Journal of Interventional Cardiology, 2009	Prospective	309	97%	4,78%	96%	NA	NA

A. Wahl	Late Results After Percutaneous Closure of Patent Foramen Ovale for Secondary Prevention of Paradoxical Embolism Using the Amplatzer PFO Occluder Without Intraprocedural Echocardiography: Effect of Device Size	JACC: Journal of the American College of Cardiology, 2008	Observational	620	100%	2.1%	91%	NA	NA
F. Wyngaert	Absence of recurrent stroke after percutaneous closure of patent foramen ovale despite residual right-to-left cardiac shunt assessed by transcranial Doppler	Archives of Cardiovascular Disease, 2008	Prospective	81	100%	NA	80%	NA	NA
B. D. Thanopoulos	Transcatheter Closure Versus Medical Therapy of Patent Foramen Ovale and Cryptogenic Stroke	Catheterization and Cardiovascular Interventions, 2006	Prospective	92	100%	29.16%	91%	NA	NA
J. M. Bijl	Percutaneous closure of patent foramen ovale	Internal Medicine Journal, 2005	Prospective	40	100%	NA	100%	NA	NA
T. Chatterjee	Interventional Closure with Amplatzer PFO Occluder of Patent Foramen Ovale in Patients with Paradoxical Cerebral Embolism	Journal of Interventional Cardiology, 2005	Prospective	55	100%	5.45%	96%	NA	NA
M. Pan	Percutaneous Treatment of Atrial Septal Aneurysm	Revista Española de Cardiología, 2005	Prospective	203	NA	3.5%	80%	NA	NA
H. W. Schuchlenz	Secondary prevention after cryptogenic cerebrovascular events in patients with patent foramen ovale	International Journal of Cardiology, 2005	Prospective	280	100%	1.11%	99%	NA	89%
A. Giardini	Comparison of Results of Percutaneous Closure of Patent Foramen Ovale for Paradoxical Embolism in Patients With Versus Without Thrombophilia	The American Journal of Cardiology, 2004	Prospective	72	100%	NA	92%	NA	NA
I. F. Purcell	Closure of Defects of the Atrial Septum in Adults Using the Amplatzer Device: 100 Consecutive Patients in a Single Center	Clinical Cardiology, 2004	Prospective	100	93%	2.13%	98%	NA	NA
H. Anzai	Incidence of Thrombus Formation on the CardioSEAL and the Amplatzer Interatrial Closure Devices	The American Journal of Cardiology, 2004	Prospective	66	100%	NA	86%	NA	NA
T. E. Hong	Transcatheter Closure of Patent Foramen Ovale Associated With Paradoxical Embolism Using the Amplatzer PFO Occluder: Initial and Intermediate-Term Results of the U.S. Multicenter Clinical Trial	Catheterization and Cardiovascular Interventions, 2003	Prospective	50	100%	NA	94%	NA	NA
K. Brockmeier	Occlusion of Interatrial Communications with the Amplatzer Device: Experience in 48 Consecutive Patients	Journal of Interventional Cardiology, 2001	Prospective	48	92%	NA	95%	NA	NA
G. Butera	Transcatheter closure of patent foramen ovale in patients with cryptogenic stroke.	Italian heart journal, 2001	Prospective	35	100%	NA	NA	NA	NA

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



* Saver JL, Carroll JD, Thaler DE, et al. Long-term outcomes of patent foramen ovale closure or medical therapy after stroke. N Engl J Med 2017; 377: 1022-32.

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