English

340027

REF V-PZ25

REF V-PZ26

Italian

REF V-PZ25

REF V-PZ26

## OC-Auto Sampling Bottle 3

## INTRODUCTION

This sampling bottle has been designed for measurement of occult blood in feces using automated feces analysis. Therefore it is ideally suited for use as a sampling bottle that can be used to easily obtain a constant amount of feces.

## PROCEDURE FOR USE

1. Twist the green cap on the sampling bottle 180 degrees to remove it.
2. The sampling probe is attached to the cap. As shown in the figure below, obtain the sample on the end of the sampling probe by scraping the probe over a broad area of the feces surface, or else insert the probe into 5 or 6 locations in the feces.
3. Return the sampling probe to the sampling bottle. Align the cap in the same direction as the bottle, then press it straight onto the bottle until a click sound is heard.
4. Shake the bottle vigorously up and down several times.

## PRECAUTIONS FOR USE

1. If the feces is hard, moisten it with water before using the sampling probe to obtain the sample.
2. If the sampling bottle (with the feces sample inside) is transported at room temperature, complete transportation within 3 days or less.
3. Do not pour out any fluid in the sampling bottle, or add water to it.
4. Do not break the aluminum seal.
5. Do not obtain feces samples during menstruation.
6. Do not use the sampling bottle by placing it in direct contact with the body.
7. Do not use the sampling bottle for any purpose other than obtaining feces samples.
8. Store the sampling bottle in a location that is out of the reach of children.
9. Perform the measurement as quickly as possible after the sample has been transported. If several days are required before the measurement is performed, store the sampling bottle at low temperature (2 - 10°C).
10. After removing the sample from the refrigerator, be sure to completely return it to room temperature before using.
11. Do not use reagents that have passed their expiration date.
12. The test sample may be contaminated with pathogens. Therefore use caution when handling.
13. After use, samples and other items must be considered to carry a risk of infection and must be treated accordingly.

## PRODUCT CODE, PRODUCT NAME, STORAGE METHOD, AND EXPIRATION DATE

Product Code	Product name	Storage method	Expiration date
V-PZ25	OC-Auto Sampling Bottle 3 with barcode	Room temp.	18 months
V-PZ26	OC-Auto Sampling Bottle 3 without barcode		

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
EC REP

AUTHORIZED REPRESENTATIVE  
WELLKANG TECH CONSULTING  
Box172 SE-15122 S-Tälje/Stockholm SWEDEN

MANUFACTURER  
EIKEN CHEMICAL CO., LTD.  
1-33-8, Hongo, Bunkyo-ku, Tokyo, 113-8408, JAPAN