330380A
330381A**OC-Control-Low**

REF V-PH57

OC-Control-High

REF V-PH58

INTRODUCTION

The amount of hemoglobin in the feces increases with illnesses that are accompanied by hemorrhagic lesions in the digestive tract, particularly in the lower digestive tract. Therefore measuring the amount of hemoglobin in the feces is an effective means of screening for early detection and treatment of colon cancer and other lower digestive tract illnesses that are accompanied by hemorrhaging.

Currently, the methods used for detection of hemoglobin include immunoassay latex agglutination, RPHA, EIA, and immunochromatography.

OC-Control-Low and OC-Control-High are freeze-dried human hemoglobin A0. They are the low- and high-concentration control samples for measurement of hemoglobin in feces. They are suitable for use as the controls for precise measurement of hemoglobin in feces.

PROCEDURE FOR USE

1. Gently remove the rubber stopper from the vial. Use a whole pipette and add precisely 1.0 mL of purified water.
2. Seal the vial with the rubber stopper, and mix by gently inverting the vial several times. Then leave for a minimum of 15 minutes.
3. Perform measurement using the same procedure as is used for a sample of unknown concentration.

REFERENCE VALUES

The attached data sheet lists the average value and reference range obtained by repeated measurements of human hemoglobin A0. The measurement values will vary somewhat depending on the skill of the measurer and other factors. Therefore it is recommended that each facility establish its own average value and reference range.

PRECAUTIONS FOR USE

1. Use the reagents on the same day they are dissolved. Store reagents frozen if it is necessary to store them after dissolution. They will remain stable for 1 month when stored at -20°C or below.
2. This product utilizes human hemoglobin A0 which is derived from human blood. The blood is tested for HBS antigens and HIV (HTLV III and LAV) antibodies, and only blood that tests negative is used. No test is done for HCV antibodies. It is impossible to completely rule out the possibility that these viruses or other pathogens may be present. Therefore, when using these reagents, handle them with the same degree of caution as patient samples.

PRODUCT CODES, PRODUCT NAMES, STORAGE METHODS, AND EXPIRATION DATES

Product Code	Product name	Storage method	Expiration date
V-PH57	OC-Control-Low	2 – 10°C	1 year
V-PH58	OC-Control-High		

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EC

REP

AUTHORIZED REPRESENTATIVE
WELLKANG TECH CONSULTING
Box172 SE-15122 S-Tälje/Stockholm SWEDEN

MANUFACTURER
EIKEN CHEMICAL CO., LTD.
33-8, Hongo 1-Chome, Bunkyo-ku, Tokyo, JAPAN

(Rev.01/2005)