

Aesculap Neurosurgery

Yasargil Aneurysm Clip System

Tray for Yasargil Aneurysm Clips



Simple, clearly organised, individual!

Yasargil Aneurysm Clip System

Practice-orientated tray for the storage, transport and sterilization of Aesculap Yasargil Aneurysm Clips.

The completed tray contains 38 storage compartments with individual identification tags.

You are free to organize the clip lay out in the tray yourself. For each clip variant there is a matching identification tag that you simply press into the tray to affix it in the desired position. This allows you to arrange the tray to correspond to your personal preferences and needs.

The identification tags contain the essential information about each clip. Their clear and uniform presentation permits fast and easy recording of the necessary clip data; helpful for classifying and reordering the clips that have been removed.

The identification tags have the same colour coding as the Yasargil Aneurysm Clips: tags for temporary clips are **yellow** and tags for permanent clips are **grey**.

Press-in tags can be individually arranged in the tray for rapid retrieval of suitable clips during surgery.



Identification tags with a uniform design containing the essential clip data and an actual size illustration of the clip.

Markings on the tray visually connect the identification tag with the corresponding clip compartment.



Transparent lid.

Secure locking system.

Length scale in mm.

Order information

Tray

Article number: **FT007P**

Size 520 x 250 x 25 mm
Autoclavable plastic
Fits 1/1 standard containers.

Identification tags

An identification tag is available for each Yasargil Aneurysm Clip. You select which clips should be in your tray and order the matching tags.

The article numbers of the tags are made up of the clip article number plus the suffix "**-TAG**".

Example:


For clip article number

FT740T



you need tag article number:

FT740T-TAG





AESCULAP®

All it takes to operate.

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen

Phone +49 7461 95-0

Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.