

Aesculap Neurosurgery

Neurosorb®


Neurosurgical Patties



In co-operation with Medthodic®

Neurosorb® –

choose from micro patties to large sheets





Size in mm

Size in inches

30x90 $1\frac{1}{6} \times 3\frac{1}{2}$

25x75 1 x 3

25x25 1 x 1

70x100 $2\frac{3}{4} \times 4$

50x90 $2 \times 3\frac{1}{2}$

Samples - Not for clinical use

Neurosorb® – your choice!

The Material

The strong pattie
very robust,
withstands rough handling.

Neurosorb 140

100% viscose

The all-round pattie
perfect balance of stability,
pliability and absorption.

Neurosorb 100

100% cotton

The green all-round
pattie
less reflection because
of green colour.

Neurosorb 100 green

100% cotton

The keyhole pattie
superior handling
preferred in micro-
neurosurgical procedures.

Neurosorb 75

100% viscose

The green keyhole
pattie
less reflection because
of green colour.

Neurosorb 75 green

100% viscose

The soft pattie
highest absorption
and excellent pliability.

Neurosorb 6

100% cotton

The delicate pattie
best pliability and high
absorption.

Neurosorb 4

100% cotton



The X-Ray Marker


no marker



dot




line



The String


no string



sewed



1/2 welded



1/1 welded



Samples – Not for clinical use

Samples - Not for clinical use

Samples – Not for clinical use

Neurosorb® – Reference Number System

1	Material
A	NS 4
B	NS 6
C	NS 75
D	NS 75 green
E	NS 100
F	NS 100 green
G	NS 140

1	X-Ray Marker
1	X-dot
2	X-line
0	no X-Mark

3-4	String
A1	sewed
B1	welded 1/1
C1	welded 1/2
Z0	no string

example

F	1	B	1	1	A	003
Neurosorb	X-dot	1/1 welded string	10 pieces	06 x 10mm		
100 green						



5	Fix
1	

6	Pieces/Set
A	10
Z	ind. Set

7-9	Size (mm)
000	04 x 06
001	04 x 10
002	06 x 06
003	06 x 10
004	06 x 15
008	10 x 10
010	10 x 20
011	10 x 25
012	10 x 30
014	10 x 50
020	15 x 15
022	15 x 25
023	15 x 30
004	15 x 40
025	15 x 55
027	15 x 75
030	20 x 20
034	20 x 50
036	20 x 60
040	25 x 25
044	25 x 75
058	30 x 90
060	05 x 150
061	15 x 150
062	06 x 80
063	06 x 40
074	50 x 90
080	70 x100



AESCULAP®

All it takes to operate.



Manufacturer according 83/42/EWG:

Medthodic

Im Süsterfeld 3
52072 Aachen / Germany
Phone +49 241 96850-0
Fax +49 241 9685056
www.medthodic.de

B | BRAUN
SHARING EXPERTISE

AESCULAP AG & CO. KG

Am Aesculap-Platz
78532 Tuttlingen / Germany
Telefon +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de