

Ordering Information

Cordis[®]
a Johnson & Johnson company

Ground breaking, Life changing™

Cardiology

cypher select[™] +

Sirolimus-eluting Coronary Stent

Make your first intervention last™



Superior Efficacy & Safety
+
Exceptional Deliverability

The CYPHER SELECT **+** Stent

Ordering Information

Crimped Profile (inch)	Nominal Stent Diameter (mm)	Nominal Stent Length (mm)					
		8	13	18	23	28	33
0.044	2.25	CRB08225	CRB13225	CRB18225	CRB23225	CRB28225	CRB33225
0.045	2.50	CRB08250	CRB13250	CRB18250	CRB23250	CRB28250	CRB33250
0.046	2.75	CRB08275	CRB13275	CRB18275	CRB23275	CRB28275	CRB33275
0.047	3.00	CRB08300	CRB13300	CRB18300	CRB23300	CRB28300	CRB33300
0.053	3.50	CRB08350	CRB13350	CRB18350	CRB23350	CRB28350	CRB33350

Technical Specifications

Stent

- Laser-cut from a 316L stainless steel tube
- Nominal wall thickness: 0.14 mm (0.0055")
- Variable strut width
- Open area is optimised around 85% for all diameters. Diameter adaptive design
- Dual markers, denoting balloon working length. The stent is positioned within the marker bands

Delivery system

- Shaft diameter: 1.9F proximal, 2.7F distal
- Nominal pressure: 11 atm (2.25 to 2.50 mm), 10 atm (2.75 to 3.50 mm)
- Rated burst pressure: 16 atm (2.25 to 3.50 mm)
- Usable length: 140 cm
- Wire compatibility: 0.36 mm (0.014")

Minimally recommended guiding catheter

- 5F (0.056" I.D.) for 2.25-3.5 mm

cypher select[™] **+**
Sirolimus-eluting Coronary Stent

Make your first intervention last[™]

Sirolimus-eluting Stent made by Cordis pursuant to a licence from Wyeth Pharmaceuticals.

Important Information: Prior to use, refer to the "instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification. The CYPHER SELECT[™] **+** Sirolimus-eluting Stent is not for sale in the USA. This publication is not intended for distribution in the USA. Europe: Tel: +32 2 352 14 11-Fax: +32 2 352 15 90 - © Johnson & Johnson Medical NV/SA • August 2006 • 2E-800-1459-1